



High Performance Refrigerator & Freezer

+4°, +5°, -20°, and -30° C

TSX Series, Cardinal Series

Service Manual

313398H06 • Revision G • 7/1/2022

Important: Read this instruction manual. Failure to follow the instructions in this manual can result in damage to the unit, injury to operating personnel, and poor equipment performance.

Caution: All internal adjustments and maintenance must be performed by qualified service personnel.

Note: Material in this manual is for informational purposes only. The contents and the product described are subject to change without notice. Thermofisher Scientific makes no representations or warranties with respect to this manual. In no event shall Thermofisher Scientific be held liable for any damages, direct or incidental, arising from, or related to the use of this manual.

Contents

| | |
|--|-----------|
| Models Covered | 5 |
| System Troubleshooting | 7 |
| Key Switch Circuit..... | 10 |
| Sensor Checking | 10 |
| FMS Unit Sensor Troubleshooting | 12 |
| 327763H01 (External Relays 1 & 2) 327325H01 (K1 – K6 and Triac #1)..... | 13 |
| External Relay 1 – Aux. Relay on Distribution Board..... | 13 |
| External Relay 2 – Light Relay on Distribution Board | 13 |
| K1 Condenser Fan Relay – High Speed..... | 13 |
| K2 Inverter Relay on 100%, Controls the Compressor Speed Through RS-232..... | 13 |
| K3 Condenser Fan Relay – Low Speed..... | 14 |
| K4 Fan Relay, Evaporator | 14 |
| Triac 1 Perimeter Heater | 14 |
| K6 Defrost / Evaporator Pan Heater Relay..... | 14 |
| 326300H01 CPU Board – DC Voltage Measurements at J12 Connector | 15 |
| 327763H01 Distribution Board (External Relays 1 & 2)..... | 16 |
| 327325H02 Relay Board (K1 – K6 and Triac 1)..... | 16 |
| 327846G02 45W Power Supply 16.5 V DC adjustable | 17 |
| 326300H01 CPU Board..... | 18 |
| TP1 “ADC0” – LM335-1 Drip Pan Sensor | 19 |
| TP2 “+5V” | 19 |
| TP3 “GND” | 20 |
| TP4 “VBAT” | 20 |
| TP5 “-5V”..... | 20 |
| TP6 “+4.5V”..... | 20 |
| TP7 “LM335-2” Lower Bottle sensor (yellow)..... | 20 |
| TP8 “+DC” | 20 |
| TP9 “LM335-3” Ambient sensor (white)..... | 20 |
| TP10 “ADC3” Cabinet Display sensor (green) | 20 |
| TP11 “ADC4” Defrost sensor (blue) | 20 |
| TP12 “ADC5” Control sensor (red)..... | 21 |
| TP14 “ADC7” AC Monitoring..... | 21 |
| TP16 “Batt Lo” | 21 |
| Position jumpers on the Control PCB: | 21 |
| Power on Reset | 21 |

| | |
|--|-----------|
| Alarm Clear / Reset..... | 21 |
| 326300H01 CPU Board Silkscreen..... | 22 |
| 326300H01 CPU Board Ladder Diagram | 23 |
| 326301H01 Display Board | 24 |
| 326301H01 Input / Output Summary..... | 24 |
| 326276H01 Display Board Silkscreen | 26 |
| 326276H01 Display Board Ladder Diagram..... | 27 |
| Cardinal Control Panel Overlays..... | 28 |
| TSX Thermo Scientific Overlays | 28 |
| Control Panel Information | 29 |
| +4° / +5° Refrigerator Alarm Functions..... | 31 |
| -20° / -30° Freezer Alarm Functions..... | 32 |
| Feature & Alarm Descriptions..... | 33 |
| Alarm Test..... | 33 |
| Audible Alarm | 33 |
| Battery Low Alarm | 34 |
| Cold Temperature Alarm..... | 34 |
| Display Low Bottle (Blood Banks only)..... | 34 |
| Door Ajar Alarm..... | 34 |
| Power Failure Alarm | 35 |
| Warm Temperature Alarm | 35 |
| Remote Alarm Contacts | 35 |
| Setpoint Security..... | 36 |
| Thermometer Gauge..... | 36 |
| Basic Functions..... | 37 |
| Programming Mode | 37 |
| Service Parameters..... | 38 |
| +4° / +5° Service Menu Functions..... | 38 |
| -20° Service Menu Functions | 40 |
| -30° Service Menu Functions | 42 |
| Firmware History | 44 |
| Parameter Defaults by Eprom Revision and Suffix | 45 |
| TSX Sensor Matrix (RoHS compliant) | 55 |
| Display Sensor and Controller Accuracy | 56 |
| Standard CPU Sensor / Harness Connections..... | 58 |
| +4°C Blood Bank Refrigerator | 58 |

| | |
|---|------------|
| +5°C Pharmacy and General-Purpose Refrigerator | 59 |
| -30°C Freezer – CFg 00 Electric Defrost | 60 |
| -30°C Freezer – CFg 01, 04 Hot Gas Defrost..... | 61 |
| -20°C ColdWall Freezer – CFg 00..... | 62 |
| -20°C 23ft ColdWall Freezer w/Bottle Sensor – CFg 02, 03, 04..... | 63 |
| -20°C 30ft ColdWall Freezer w/Bottle Sensor – CFg 01, 02, 03, 04..... | 64 |
| Sensor Locations Inside of Various TSX Units..... | 65 |
| Inverter Information | 69 |
| Embraco High Performance (HP) | 69 |
| Embraco High Performance (HP) LED Codes | 71 |
| Maia Inverter (Started with some Suffix 17 Models)..... | 72 |
| Maia LED Codes | 73 |
| Condenser Fan..... | 74 |
| Electric Condensate Heater 328840H0X..... | 75 |
| Chromatography with Safety Outlet..... | 76 |
| HC Sensor PPM Sensativity Graph..... | 79 |
| Safety Circuit Wiring Diagram (328790H01) | 80 |
| Blood Bank Anti-Freeze Protection | 81 |
| Dixell Operating Manual..... | 83 |
| Anti-Tip Bracket Instructions..... | 85 |
| 4-20mA Temperature Transmitter Information | 86 |
| 4-20mA Option Numbers | 89 |
| 4-20mA Temperature Transmitter Wiring Diagram (327987H01)..... | 90 |
| Wiring Diagram +4°C Refrigerator Suffix 17 / 18 328059H05..... | 91 |
| Wiring Diagram +5°C Refrigerator Suffix 17 / 18 328059H08..... | 92 |
| Wiring Diagram -20°C Coldwall Suffix 17 – 20 Freezer 328059H06 | 93 |
| Wiring Diagram -30°C Suffix 17 / 18 Freezer 328059H07 | 94 |
| Charge Information..... | 95 |
| Pressure and Amp Draw | 96 |
| Refrigeration Diagram +4°C / +5°C Suffix 17, 18 327990H04..... | 99 |
| Refrigeration Diagram -20°C ColdWall Suffix 17 - 20 327990H06 | 100 |
| Refrigeration Diagram -30°C Suffix 17 / 18 327990H05 | 101 |
| Touch Up Paint Information..... | 102 |
| Technical Service Bulletins of Note..... | 103 |
| Available Service Kits..... | 104 |

Models Covered

TSX Models covered by this Service Manual

| | | | | |
|-----------|-----------|-----------|-----------|-----------|
| TSX1204BA | TSX2304BZ | TSX3004BY | TSX5004BV | TSX4505CD |
| TSX1204BD | TSX2305CA | TSX3004BZ | TSX5004BY | TSX4505CV |
| TSX1204BV | TSX2305CD | TSX3005CA | TSX5004BZ | TSX4505CY |
| TSX1204BY | TSX2305CV | TSX3005CD | TSX5005CA | TSX4505CZ |
| TSX1204BZ | TSX2305CY | TSX3005CV | TSX5005CD | TSX4505GA |
| TSX1205GA | TSX2305CZ | TSX3005CY | TSX5005CV | TSX4505GD |
| TSX1205GD | TSX2305GA | TSX3005CZ | TSX5005CY | TSX4505GV |
| TSX1205GV | TSX2305GD | TSX3005GA | TSX5005CZ | TSX4505GY |
| TSX1205GY | TSX2305GV | TSX3005GD | TSX5005GA | TSX4505GZ |
| TSX1205GZ | TSX2305GY | TSX3005GV | TSX5005GD | TSX2320HA |
| TSX1205PA | TSX2305GZ | TSX3005GY | TSX5005GV | TXS2320EA |
| TSX1205PD | TSX2305PA | TSX3005GZ | TSX5005GY | TXS2320ED |
| TSX1205PV | TSX2305PD | TSX3005PA | TSX5005GZ | TXS2320EV |
| TSX1205PY | TSX2305PV | TSX3005PD | TSX5005PA | TXS2320EY |
| TSX1205PZ | TSX2305PY | TSX3005PV | TSX5005PD | TXS2320EZ |
| TSX1205SA | TSX2305PZ | TSX3005PY | TSX5005PV | TXS2320FA |
| TSX1205SD | TSX2305SA | TSX3005PZ | TSX5005PY | TXS2320FD |
| TSX1205SV | TSX2305SD | TSX3005SA | TSX5005PZ | TXS2320FV |
| TSX1205SY | TSX2305SV | TSX3005SD | TSX5005SA | TXS2320FY |
| TSX1205SZ | TSX2305SY | TSX3005SV | TSX5005SD | TXS2320FZ |
| TSX1230FA | TSX2305SZ | TSX3005SY | TSX5005SV | TXS3020EA |
| TSX1230FD | TSX2330FA | TSX3005SZ | TSX5005SY | TXS3020ED |
| TSX1230FV | TSX2330FD | TSX3030FA | TSX5005SZ | TXS3020EV |
| TSX1230FY | TSX2330FV | TSX3030FD | TSX5030FA | TXS3020EY |
| TSX1230FZ | TSX2330FY | TSX3030FV | TSX5030FD | TXS3020EZ |
| TSX1230LA | TSX2330FZ | TSX3030FY | TSX5030FV | TXS3020FA |
| TSX1230LD | TSX2330LA | TSX3030FZ | TSX5030FY | TXS3020FD |
| TSX1230LV | TSX2330LD | TSX3030LA | TSX5030FZ | TXS3020FV |
| TSX1230LY | TSX2330LV | TSX3030LD | TSX5030LA | TXS3020FY |
| TSX1230LZ | TSX2330LY | TSX3030LV | TSX5030LD | TXS3020FZ |
| TSX2304BA | TSX2330LZ | TSX3030LY | TSX5030LV | |
| TSX2304BD | TSX3004BA | TSX3030LZ | TSX5030LY | |
| TSX2304BV | TSX3004BD | TSX5004BA | TSX5030LZ | |
| TSX2304BY | TSX3004BV | TSX5004BD | TSX4505CA | |

CAX Models covered by this Service Manual

| | | | |
|----------|----------|----------|----------|
| CH1204BA | CH2305PA | CH3005PA | CH5004PA |
| CH1204BD | CH2305PD | CH3005PD | CH5004PD |
| CH1205GA | CH2305SA | CH3005SA | CH5004SA |
| CH1205GD | CH2305SD | CH3005SD | CH5004SD |
| CH1205PA | CH2330FA | CH3030FA | CH5030FA |
| CH1205PD | CH2330FD | CH3030FD | CH5030FD |
| CH1205SA | CH2330LA | CH3030LA | CH5030LA |
| CH1205SD | CH2330LD | CH3030LD | CH5030LD |
| CH1230FA | CH2320EA | CH3020EA | CH4505CA |
| CH1230FD | CH2320ED | CH3020ED | CH4505CD |
| CH1230LA | CH2320FA | CH3020FA | CH4505GA |
| CH1230LD | CH2320FD | CH3020FD | CH4505GD |
| CH2304BA | CH3004BA | CH5004BA | |
| CH2304BD | CH3004BD | CH5004BD | |
| CH2305CA | CH3005CA | CH5005CA | |
| CH2305CD | CH3005CD | CH5005CD | |
| CH2305GA | CH3005GA | CH5004GA | |
| CH2305GD | CH3005GD | CH5004GD | |

System Troubleshooting

| | |
|---|---|
| Unit is Completely Non-Functional | <ul style="list-style-type: none"> • Check the circuit breaker on the back of the unit is in the ON position • Check power cord and outlet for proper line voltage • Check that the key switch position is in the On or Alarm positions • Check for shorted sensors – see Sensor Checking following this chart • Check for 16.5v DC at J12 pin 1 of CPU board – if not present, replace 45W Power Supply board • Check for PS15R DC return signal through the key switch – see “Key Switch Circuit” • Check for 5v DC at J12 pin 2 of CPU board – if not present, replace CPU board • Check for proper orientation of the EPROM (if power has been applied while the chip was inserted 180° out, replace the EPROM) • Replace CPU / Power Relay Board with appropriate kit. • Replace 45W Power Supply board • Replace Display board • Replace Distribution board |
| Unit does not operate; thermometer bulb is the only thing displayed | <ul style="list-style-type: none"> • System cannot acquire a good read from the cabinet sensor; cables are too close to power lines, fans, or other interference sources, or a strong RF source is nearby. Reroute sensor cables. • Check sensors for open / short – see Sensor Checking following this chart • Replace CPU / Power Relay Board with appropriate kit. |
| Display reads all 8's | <ul style="list-style-type: none"> • Check connections between the CPU board and the display board • Replace RJ45 6 Pin cable • Replace display board |
| Unit runs 100% and is warmer than set point | <ul style="list-style-type: none"> • Check refrigeration system if unit is not reaching set point • Check for a bad door gasket or port hole air leaks |
| Unit runs 100% and is colder than set point | <ul style="list-style-type: none"> • Check system calibration by comparing setpoint, display, and measured chamber temperature • Check for loose or open cabinet sensor – see Sensor Checking following this chart • Replace CPU / Power Relay Board with appropriate kit. • Replace 45W Power Supply board • Replace Distribution board |
| Compressor will not start even though cabinet temperature is warmer than setpoint | <ul style="list-style-type: none"> • Unit could be in short cycle period, wait 2 minutes for normal cycle and 7 minutes after a power failure to verify • If +4° / +5°C unit, check control probe reading in Service Mode – see 4° Operation • Verify HPCO jumper is installed on the Relay board if this is not a 3030, or 5030 freezer |

| | |
|--|---|
| “Err” Displayed | <ul style="list-style-type: none"> • “Err” indicates the CPU can not read an ADC channel – applies to sensors and line voltage reading – see Sensor Checking following this chart • If occurs during normal run mode, replace the display sensor (green) • If occurs while reading the lower bottle (blood banks only), replace lower bottle sensor (yellow) • If occurs while reading the defrost sensor in service mode, then replace the defrost sensor (blue) • If occurs while reading the control sensor in service mode, replace the control sensor (red) • If occurs while reading the ambient sensor in service mode, replace the ambient sensor (white) • If occurs while reading the drip pan sensor in service mode, calibrate R22, replace the drip pan sensor (black) • Replace CPU / Power Relay Board with appropriate kit. |
| Display will not respond to keypad inputs | <ul style="list-style-type: none"> • Loose or incorrect connection between front panel overlay and display board – harness off by one pin • Check J3 on display board – pin 4 (bottom) is ground; other three should be at +5v DC until a button is pressed; pressing “down arrow” should make pin 3 go to 0v; pressing center button should make pin 2 go to 0v, pressing “up arrow” should make pin 1 go to 0v. If any one does not work, replace the overlay • Replace the RJ45 6 Pin and 10 Pin cables • Replace the Display board • Replace CPU / Power Relay Board with appropriate kit. |
| Button on overlay is acting as if it is stuck | <ul style="list-style-type: none"> • Insulate the overlay cable strip from the front panel |
| Unit always times out of defrost | <ul style="list-style-type: none"> • Defective RTD defrost sensor, check valued in Service Menu, replace if “Err” or reading over 40°C • Replace CPU / Power Relay Board with appropriate kit. |
| Evaporator fan is not turning (with door closed and unit not in defrost) | <ul style="list-style-type: none"> • Check door switch operation; Normally Closed – held open by the door. Door harness at J9 on CPU board can be removed to simulate a closed door, if the fan comes on, check the door switch wiring and door alignment • Check defrost sensor reading in Service Mode, must be below –1°C for freezers • Check for line voltage at the fan motor(s) • Replace CPU / Power Relay Board with appropriate kit. • Replace the Distribution board |
| Unit does not go into Power Fail mode during loss of line voltage | <ul style="list-style-type: none"> • Weak battery – replace the battery. Limited shelf life of 6 months on 3.3 amp hour batteries. • Replace CPU / Power Relay Board with appropriate kit. • Replace 45W Power Supply Board |
| Unit goes into Power Fail mode, but there is no loss of line voltage | <ul style="list-style-type: none"> • Replace CPU / Power Relay Board with appropriate kit. |

| | |
|---|---|
| Interior lights do not illuminate with the front panel light switch on | <ul style="list-style-type: none"> • Check the light switch connections at the front panel • Check / replace the LED power supply • Check / replace the LED lights • Replace CPU / Power Relay Board with appropriate kit. |
| Interior lights work with the front panel switch, but do not illuminate with the door open | <ul style="list-style-type: none"> • Short CPU board J9 pins 1 & 2 to simulate an open door – if the lights illuminate, check door switch(s) and harness connections • Replace CPU / Power Relay Board with appropriate kit. |
| Interior lights stay on even with the door(s) closed and the front panel switch in the off position | <ul style="list-style-type: none"> • Check door alignment , may be sagging causing the switch to stay active • Check the light switch for a short • Replace CPU / Power Relay Board with appropriate kit. |
| Showing E01 on Display - Invalid Mode in the Code | <ul style="list-style-type: none"> • New set of boards installed and three parameters in the Service Menu need to be setup to identify what kind of unit you are working with before the software will run – update Type, Config, & Cab Size, exit service, then wait for the unit to reboot itself to accept the changes. • <i>Do not go onto the following steps before the previous setup is confirmed by going back into the Service Menu</i> • Replace the RJ45 cable between the CPU and Relay boards • Replace the RJ45 cables between the CPU and Display boards • Replace CPU / Power Relay Board with appropriate kit. |
| Showing E02 on Display - Control sensor failure | <ul style="list-style-type: none"> • Replace control sensor (red) • Replace CPU / Power Relay Board with appropriate kit. |
| Showing E03 on Display - Defrost sensor failure | <ul style="list-style-type: none"> • Replace defrost sensor (blue) • Replace CPU / Power Relay Board with appropriate kit. |
| Showing E05 on Display - Ambient sensor failure | <ul style="list-style-type: none"> • Replace ambient sensor (white) • Replace CPU / Power Relay Board with appropriate kit. |
| Showing E07 on Display - Compressor running more than 2 hours | <ul style="list-style-type: none"> • Check the PEP for air blockage • Defrost the coldwall unit |
| Showing “----“ on Display- Lost Communications | <ul style="list-style-type: none"> • Check / Replace RJ45 cable between CPU and Distribution boards • Replace CPU / Power Relay Board with appropriate kit. • Replace Distribution board |

Key Switch Circuit

With the unit plugged in, there should be 16.5v DC (signal named - **PS15**) at CPU **J12-8** provided by the 45w Power Supp board through Distribution Board **J9-3**.

Once the key switch is turned on, you should get 16.5v DC (**PS15R**) back on **J12-9**, which jumpers back to **J12-10** to power up the CPU board.

If you are not getting **PS15R** back, check along the path for an open circuit.

PS15 signal leaves the Distribution board **J9-3** to the CPU board **J12-8**; then out of the CPU board **J7-1** to the display board **J2-1**; then out of the display board **J4-4** to key switch terminal **#5**, then out key switch terminals **#7** or **#8** (**depends if in the On or Alarm position**) as **PS15R** back to display board **J4-5** and **J4-6**. **PS15R** leaves the display board **J2-10** to CPU board **J7-10**; then out of the CPU board **J12-9** to CPU board **J12-10**.

NOTE – the CPU and display boards just pass the signals **PS15** and **PS15R** through them, so an open circuit would most likely be at the connection point between the boards and harnesses, in one of the harnesses, or in the key switch.

Sensor Checking

Refer to the figures below for pin outs and the following RTD Temp/Ohms Chart for 1000 ohm sensors. 1000 ohm sensors are connected with redundant wires for reliability. **NOTE** - An easy way to read the pins on the sensors is to purchase a Cat 5 network jack from any electronics or home improvement store that handles networking products. The sensors use an 8 wire network cable, so the jack will extend the connections for easier access with meter leads.

309159H09 (Green 72") +4°, +5°, -20° / -30° Cabinet sensor – pins 1,5,6 to 2,3,4 should be the RTD value, pins 7 to 8 should be a fixed resistor of 235 ohms +/- 20 ohms.

309159H10 (green 120") +4°, +5°, -30° Cabinet sensor – pins 1,5,6 to 2,3,4 should be the RTD value, pins 7 to 8 should be a fixed resistor of 235 ohms +/- 20 ohms.

326737H05 (blue) +4°, +5°, -30° defrost sensor – pins 1,5,6 to 2,3,4 should be the RTD value, pins 4 to 8 should be an open circuit.

326737H03 (red) +4, -20°, -30° Control sensors – pin 1,5,6 to 2,3,4 should be the RTD value, pin 7 & 8 should be open.




326737H04 (red) +4, -20°, -30° Control sensors – pin 1,5,6 to 2,3,4 should be the RTD value, pin 7 & 8 should be open.

326739H01 (white) LM335 Ambient sensor can only be checked in an active circuit.

312096H06 (yellow) +4° lower bottle sensor can only be checked in an active circuit.

326739H03 (black) -30° drip pan sensor can only be checked in an active circuit – TP1

85089G10 (gray) 4-20mA and chart recorder sensors - wire 3 & 4 should be the **100** ohm RTD value. (you can still use the chart, just divide the 1000 ohm resistance value by 10 to give the approximate 100 ohm RTD value)





NOTE: Sensors with any connections between the RTD pins and pins with other devices (fixed resistor, shield ground or TE devices) are considered defective sensors.

RTD Temperature in °C Vs. Resistance in Ohms

| Temperature | Resistance | Temperature | Resistance | Temperature | Resistance |
|-------------|------------|-------------|------------|-------------|------------|
| -50 | 803.07 | -21 | 917.66 | 8 | 1031.22 |
| -49 | 807.03 | -20 | 921.59 | 9 | 1035.12 |
| -48 | 811.00 | -19 | 925.53 | 10 | 1039.02 |
| -47 | 814.97 | -18 | 929.46 | 11 | 1042.92 |
| -46 | 818.93 | -17 | 933.38 | 12 | 1046.81 |
| -45 | 822.90 | -16 | 937.31 | 13 | 1050.71 |
| -44 | 826.86 | -15 | 941.24 | 14 | 1054.60 |
| -43 | 830.82 | -14 | 945.16 | 15 | 1058.49 |
| -42 | 834.78 | -13 | 949.09 | 16 | 1062.38 |
| -41 | 838.74 | -12 | 953.01 | 17 | 1066.27 |
| -40 | 842.70 | -11 | 956.93 | 18 | 1070.16 |
| -39 | 846.66 | -10 | 960.85 | 19 | 1074.04 |
| -38 | 851.61 | -9 | 964.77 | 20 | 1077.93 |
| -37 | 854.57 | -8 | 968.69 | 21 | 1081.81 |
| -36 | 858.52 | -7 | 972.61 | 22 | 1085.70 |
| -35 | 862.47 | -6 | 976.52 | 23 | 1089.58 |
| -34 | 866.42 | -5 | 980.44 | 24 | 1093.46 |
| -33 | 870.37 | -4 | 984.35 | 25 | 1097.34 |
| -32 | 874.32 | -3 | 988.26 | 26 | 1101.22 |
| -31 | 878.27 | -2 | 992.18 | 27 | 1105.10 |
| -30 | 882.21 | -1 | 996.09 | 28 | 1108.97 |
| -29 | 886.16 | 0 | 1000.00 | 29 | 1112.85 |
| -28 | 890.10 | 1 | 1003.90 | 30 | 1116.73 |
| -27 | 894.04 | 2 | 1007.81 | 31 | 1120.61 |
| -26 | 897.98 | 3 | 1011.71 | 32 | 1124.49 |
| -25 | 901.92 | 4 | 1015.62 | 33 | 1128.37 |
| -24 | 905.86 | 5 | 1019.52 | 34 | 1132.24 |
| -23 | 909.79 | 6 | 1023.42 | 35 | 1136.12 |
| -22 | 913.73 | 7 | 1027.32 | | |

FMS Unit Sensor Troubleshooting

FMS (Flammable Storage) Units per UL requirements the sensors must be “Intrinsically Safe”. This introduced the need for Barriers in the sensor circuits. These barriers effect the sensor offsets and troubleshooting.



Wires and barriers are labeled 1 thru 4

Note: To do sensor checking the sensor needs to be disconnected from the intrinsic Barrier.

327763H01 (External Relays 1 & 2) 327325H01 (K1 – K6 and Triac #1)

External Relay 1 – Aux. Relay on Distribution Board

[4°, 5°, -20°, -30°] Under normal operation, this relay will provide power to the drain tube heater(if used).

External Relay 1 should be on whenever there is power to the unit, the key switch is on, and the software has started.

External Relay 2 – Light Relay on Distribution Board

[4°, 5°] On if door is open. Off if door is closed. The front panel light switch forces this relay ON regardless of door position.

[-20°, -30°] **Never turns on; no interior lights.**

K1 Condenser Fan Relay – High Speed

[4°] Turns On if control probe (evaporator airflow) reads a temperature \geq [Set point + 'UCL'] value and minimum off time requirement has been met. Turns on 10 seconds before the compressor. Fan control may revert to Low speed (K3 on at the same time) if that is all that is required to maintain cabinet temp.

Turns OFF when control probe reads a temperature $<$ [Set point – 'LCL'] and the minimum on time requirement has been met. Minimum on time is approximately 3 minutes and is “hard coded” and is not adjustable. Minimum off time is approximately 6 minutes and is “hard coded” and is not adjustable to prevent short cycling. Turns off 10 seconds after the compressor has stopped.

[-20°, -30°] Turns On if red control probe reads a temperature \geq [Set point + UCL] and minimum off time requirement has been met. Turns on 10 seconds before the compressor. Fan control may revert to Low speed (K3 on at the same time) if that is all that is required to maintain cabinet temp.

Turns Off if red control probe reads a temperature $<$ [Set point – LCL] and minimum on time requirement has been met. Minimum on time is approximately 3 minutes. Minimum off time is approximately 4-6 minutes. Minimum values are “hard coded” and are not adjustable. Turns off 10 seconds after the compressor has stopped.

K2 Inverter Relay on 100%, Controls the Compressor Speed Through RS-232

[4°, 5°] Compressor turns On if control probe (evaporator airflow) reads a temperature \geq [Set point + 'UCL'] value and minimum off time requirement has been met. Turns OFF when control probe reads a temperature $<$ [Set point – 'LCL'] and the minimum on time requirement has been met. Minimum on time is approximately 2 minutes and is “hard coded” and is not adjustable. Minimum off time is approximately 6 minutes and is “hard coded” and is not adjustable to prevent short cycling.

[-20°, -30°] Compressor turn on if red control probe reads a temperature \geq [Set point + UCL] and minimum off time requirement has been met. Turns Off if red control probe reads a temperature $<$ [Set point – LCL] and minimum on time requirement has been met. Minimum on time is approximately 2 minutes. Minimum off time is approximately 4-6 minutes. Minimum values are “hard coded” and are not adjustable.

K3 Condenser Fan Relay – Low Speed

[4°, 5°] Turns On (along with K1) if control probe (evaporator airflow) reads a temperature \geq [Set point + ‘UCL’] value and minimum off time requirement has been met. Turns on 10 seconds before the compressor. Fan control may revert to just K1 on - High speed if that is required to maintain cabinet temp.

Turns OFF when control probe reads a temperature $<$ [Set point – ‘LCL’] and the minimum on time requirement has been met. Minimum on time is approximately 2 minutes and is “hard coded” and is not adjustable. Minimum off time is approximately 6 minutes and is “hard coded” and is not adjustable to prevent short cycling. Turns off 10 seconds after the compressor has stopped.

[-20°, -30°] On (along with K1) if red control probe reads a temperature \geq [Set point + UCL] and minimum off time requirement has been met. Turns on 10 seconds before the compressor. Fan control may revert to K1 High speed if that is required to maintain cabinet temp. Turns Off if red control probe reads a temperature $<$ [Set point – LCL] and minimum on time requirement has been met. Minimum on time is approximately 2 minutes.

Minimum off time is approximately 4-6 minutes. Minimum values are “hard coded” and are not adjustable. Turns off 10 seconds after the compressor has stopped.

K4 Fan Relay, Evaporator

[4°, 5°] If there is power to the unit, the key switch is on and the software has started, this relay should always be on. Door openings will turn the fan off, but only for 6 minutes. Then the fan(s) turn back on regardless of door position.

[-20°] **Never turns on; no evaporator fan on this model.**

[-30°] Turns on once defrost probe reads a temperature $<$ -20°C. Turns off during a defrost cycle; back on after the defrost cycle and when defrost probe reads a temperature $<$ -20°C. Turns off if the door is open, but only for 6 minutes, then the fan(s) are forced back on regardless of door positions. Turns on once the door is closed, if the defrost probe reads a temperature $<$ -20°C.

Note: The Firmware has built in hidden timers that operate in the following way:

They will start the fans at 60 minutes during pulldown or 10 minutes during normal cycling if there is an issue with the defrost probe reaching $<$ -20°C. In the case of defrost probe failure the fans will start without waiting.

Triac 1 Perimeter Heater

[4°, 5, -20°, -30°] The perimeter heater is now controlled with pulse with modulation (percent on time). The heater percentage “on time” is preset based on the program default table found later in the service manual. The setting is adjustable through the Service Menu with a range of 20% to 100% on time.

K6 Defrost / Evaporator Pan Heater Relay

[4°, 5°, -20°] Never turns on; no heater on these models.

[-30°] Turns on approximately six hours after the last defrost cycle was terminated (six hour count is initiated at the termination of a defrost cycle or if unit is cycled off with the key). Turns off after 20 minutes or when defrost probe reads a temperature $>$ +1°C.

NOTE: The Defrost parameters are now settable in the Service Menu if **Adaptive Defrost** is disabled, so those settings will alter Interval, Time out and Temp out.

326300H01 CPU Board – DC Voltage Measurements at J12 Connector

| J12 Pin Number | Signal Name | Wire Color | Power Applied – Key Switch Off | Power Applied – Key Switch ON / Alarm | Notes |
|----------------|-------------|-------------|--------------------------------|---------------------------------------|---|
| 1 | DC In | Red | 16.5v DC | 16.5v DC | From 45W Power Supply |
| 2 | +5v DC | Green | 0v DC | 5.0v DC | Supplied by CPU Board |
| 3 | NC | | NA | NA | No Connection |
| 4 | NC | | NA | NA | No Connection |
| 5 | NC | | NA | NA | No Connection |
| 6 | NC | | NA | NA | No Connection |
| 7 | NC | | NA | NA | No Connection |
| 8 | PS15 | Red | 16.5v DC | 16.5v DC | Power Supply to Key Switch |
| 9 | PS15R | Red / White | 0v DC | 16.2v DC | Jumper to Pin 10 (Return from Key Switch) |
| 10 | +12v DC | Red / White | 0v DC | 16.2v DC | Jumper to Pin 9 (Powers Up 12v Signal) |
| 11 | Ground | Brown | 0v DC | 0v DC | Should equal frame ground |
| 12 | NC | | NA | NA | No Connection |
| 13 | NC | | NA | NA | No Connection |

NOTE: Any time a unit is serviced, it should be made a practice to check that the nubbin is knocked off the upper left of each relay. This can be accomplished quickly & simply with a flat screwdriver. This provides ventilation for the relay, reducing hydrocarbon buildup, and ultimately increases the life of the relay.


|----- Relays and Triac that reside on the Relay Bd -----| |---Distribution BD ---|

| Model | K1 (Cond. Fan) | K2 (Inverter) | K3 (Cond. Fan Low Speed) | K4 (Evap Fan) | Triac -1 (Perimeter heater) | K6 (Defrost Heater) | Ext. Rly - 1 (Aux. ON) | Ext Rly - 2 (Light) |
|-------|-------------------|------------------|--------------------------------|------------------|-----------------------------------|---------------------------|---------------------------|------------------------|
| +4C | √ | √ | √ | √ | √ | X | √ | √ |
| -20C | √ | √ | √ | X | √ | X | √ | X |
| -30C | √ | √ | √ | √ | √ | √ | √ | X |

327763H01 Distribution Board (External Relays 1 & 2)



327325H02 Relay Board (K1 – K6 and Triac 1)



327846G02 45W Power Supply 16.5 V DC adjustable

Accepts input line voltage range of 90v to 264v AC 50/60 Hz.




Basically a 15v DC supply that should be preadjusted to 16.5v DC.

Use 327846G02 S service kit to replace the former H01 boards. Kit contains the new H02 smaller supply, adapter plate to match the H01 mounting, and two harnesses.

Note: There are no static sensitive components on this board, so replacement may arrive without the typical anti static bag.

326300H01 CPU Board



Test Points (as labeled on the silkscreen)

| | | |
|--------------|---------------|-------------------|
| TP1 → "ADCO" | TP6 → "+4.5V" | TP11 → "ADC4" |
| TP2 → "+5V" | TP7 → "ADC1" | TP12 → "ADC5" |
| TP3 → "GND" | TP8 → "+DC" | TP13 → "ADC6" |
| TP4 → "VBAT" | TP9 → "ADC2" | TP14 → "ADC7" |
| TP5 → "-5V" | TP10 → "ADC3" | TP15 → "-1.6 Adj" |
| | | TP15 → "Batt Lo" |

General Overview

| Test Point | Generic Description | H01 Application |
|------------|---|-------------------------|
| TP1 "ADCO" | Voltage @ input to the A/D converter (channel 0). aka LM335-1 | Drip Pan Sensor (Black) |
| TP2 "+5V" | 5 VDC, regulator output. | 5VDC |
| TP3 "GND" | Ground plane. | GND |
| TP4 "VBAT" | Positive terminal of the battery. Charging voltage under normal operation; battery voltage in backup. | VBAT |
| TP5 "-5V" | -5 VDC, regulator output. | Not Used |
| TP6 "4.5V" | +4.5VDC, regulator output. | +4.5VDC |

| Test Point | Generic Description | H01 Application |
|--------------------|---|-----------------------------------|
| TP7 “ADC1” | Voltage @ input to the A/D converter (channel 1). aka LM335-2 | Lower Bottle Sensor |
| TP8 “+DC” | 12-15 Volt DC Input to the 5v DC regulator | Input to the 5v DC Regulator |
| TP9 “ADC2” | Voltage @ input to the A/D converter (channel 2). aka LM335-3 | Ambient Sensor (White) |
| TP10 “ADC3” | Voltage @ input to the A/D converter (channel 3). aka RTD-1 | Display Sensor (Green) |
| TP11 “ADC4” | Voltage @ input to the A/D converter (channel 4). aka RTD-2 | Defrost Sensor (Blue) |
| TP12 “ADC5” | Voltage @ input to the A/D converter (channel 5). aka RTD-3 | Control Sensor (Red) |
| TP13 “ADC6” | Voltage @ input to the A/D converter (channel 6). | Not Used (Was Line Voltage) |
| TP14 “ADC7” | Voltage @ input to the A/D converter (channel 7). | Monitor AC Input for Power Fail |
| TP15 “-1.6 Adj.” | Offset voltage adjustment specific to ULT T/C circuitry. | Not Used |
| TP16 “Battery Low” | Comparator circuit for batter level Low (~0.095v) = OK, High (~5v) = low battery | Used to Monitor the Battery Level |

Note:

1. When observing the voltage on any A/D channel inputs, there will be a voltage ‘pulse’ that occurs. This is normal and should not cause alarm.
2. If **any channel** is found to be outside the limits of normal operation, 0 to 4.5 VDC, the reading of all channels will become suspect. For instance, if channel 7 is found to be 5.4 VDC, the readings on channels 0-6 may be corrupt. Check the sensor for the corresponding channel, if sensor is not bad, replace the CPU board.

TP1 “ADC0” – LM335-1 Drip Pan Sensor

Black sensor used to terminate pan heat after the defrost cycle is complete. Basically controls how long drip time is by temping out before the max of 5 minutes is up. Drip pan is also called the evaporator cover.

TP2 “+5V”

CPU board converts +DC (16.5v) voltage from Power Supply board into 5v DC to drive all of the chips on the CPU board. 5v not present will cause a “dead” system. Check between pins 10 & 11 of J12, if 12-15v DC is present, replace CPU board, otherwise, check harness / Power Supply board.

TP3 “GND”

Frame ground.

TP4 “VBAT”

When operating correctly, 14 – 16ish VDC is typically present when loaded with a battery. 17 – 19ish VDC will be present when unloaded.

TP5 “-5V”

For the -86°C TC circuit – NOT use in this application.

TP6 “+4.5V”

Reference voltage for the sensors. Incorrect voltage reading will effect RTD sensor readings. Since 4.5v is generated on CPU board – check for shorted sensors, then replace CPU board.

TP7 “LM335-2” Lower Bottle sensor (yellow)

Lower Bottle sensor (was 1000 ohm RTD), now requires an active circuit due to being a LM335 IC

TP8 “+DC”

DC voltage from the Power Supply board, will read in the 12 –15v range based on actual line voltage. Used to generate 5v DC and 4.5v reference. If the voltage is higher than 18v or less than 11v, check the voltage program jumpers on the Power Supply board. If the jumper is not set to match the actual line voltage, reposition it. Otherwise, replace the Power Supply board.

TP9 “LM335-3” Ambient sensor (white)

Ambient sensor changed from previous defrost role, but still requires an active circuit due to being a LM335 IC

TP10 “ADC3” Cabinet Display sensor (green)

Upper Cabinet 1000 ohm RTD sensor located in solution bottle (except minus 20 cold wall). Used for control / display on the –20° and –30° units. Used just for display on the +4° units. Should read in the 0.082v (freezers) to 0.095v (refrigerators) range. Open circuit will read in the 4.5v area. Short circuit will read in the 0v area.

TP11 “ADC4” Defrost sensor (blue)

Defrost sensor in the evaporator 1000 ohm RTD sensor. Open circuit will read in the 4.5v area. Short circuit will read in the 0v area.

TP12 “ADC5” Control sensor (red)

Control 1000 ohm RTD sensor located in the air in front of the defrost sensor inside the evaporator in refrigerators. Should read in the 0.095v area. Open circuit will read in the 4.5v area. Short circuit will read in the 0v area. On the side wall in -20 / -30° Freezers.

TP14 “ADC7” AC Monitoring

Used to monitor the AC power for power fail conditions. Should read in the 1.8v area.

TP16 “Batt Lo”

Output of the battery level comparator circuit attached to the VBatt signal. Normal operation while the battery is charging or at a voltage over 9.6V, the output will be less than 1v. When the battery is tested under load and the voltage drops below 9.6v, the output will change to greater than 4v.

Position jumpers on the Control PCB:

- **JP1** - not on NGC3 board set (was 50/60Hz setting)
- **JP2** - connect pins 1 & 2 for 256k PROM chip
- **JP3** - connect pins 1 & 2 for Telemode (display will be blank unless the jumper is placed in this configuration)
- **JP20** - connect pins 2 & 3 for 5891 over-ride not active

Power on Reset

Press Up and Down arrows together during power up sequence to reset the on-board NVRAM to factory default programming and service settings.

CAUTION


Power on Resets are necessary when changing EEPROMS, and sometimes used to correct service issues, but a Power on Reset **will overwrite** all settings customized to the unit (cut on, differential, etc.). **So, write them down first!**

Note:


Changes to **Type, Config, Cu ft** parameters causes the software to reboot or POR. Reboot may take up to **3 minutes**.

Alarm Clear / Reset


Hold up and down arrows together during normal operation to reset alarms, or “Has Been” Alarms, plus cold and warm temperature excursions (max/min) stored values. Alarm Ringback if it has not been corrected is 10 minutes.



326300H01 CPU Board Ladder Diagram



326301H01 Display Board



326301H01 Input / Output Summary

J2 – General I/O with CPU board

- Pin 1 – “PS15” DC voltage supplied out to key switch. Refer to **J4** for troubleshooting.
- Pin 2 – “+12V” DC voltage supplied to audio transducer (AL1), 5 volt regulator (U2), and the VFD in ULT applications. The voltage level at “+12V” should be one diode drop down from “+DC” (CPU board and Power Supply board). Typically, these will be 10.5ish VDC and 11.3ish VDC, respectively. Can also be check at ‘+’ side of AL1 or U2 pin 1 (toward J5 & J2).
- Pin 3 – “SW-DOWN” Keypad status passed to CPU board. See **J3** for troubleshooting.
- Pin 4 – “SW-SELECT” Keypad status passed to CPU board. See **J3** for troubleshooting.
- Pin 5 – “SW-UP” Keypad status passed to CPU board. See **J3** for troubleshooting.
- Pin 6 – Ground
- Pin 7 – “BEEPER” Signal from CPU board to activate an audible alarm. +5 VDC when active, measure at the anode (+) of D1 (pin away from audio transducer - AL1). Anode is also the side of the diode without a bar.
- Pin 8 – “ALARM-EN” Output supplied to CPU board indicating that the key switch is in position 3. +5 VDC when active.
- Pin 9 – “BATTERY-EN” Output supplied to battery charging circuit. Refer to **J4 pin 2 & 3** for troubleshooting.
- Pin 10 – “PS15R” Secondary side of key switch. Refer to **J4** for troubleshooting.

J3 – Switch Input from Keypad

Pressing a keypad simply shorts it to ground. For instance, pressing the UP pad will cause pin 1 to short to pin 4. Pin 1 is toward the top of the board or R2/D1. **NOTE Pin 1 on a connector can frequently be identified by a square pad on the PCB.**

Pin 1 – “SW-UP” Input from UP keypad. Grounding this pin will simulate the ‘up’ pad being pressed.

Pin 2 – “SW-SELECT” Input from Mode/Mute/Select keypad. Grounding this pin will simulate the ‘Mode/Mute/Select’ pad being pressed.

Pin 3 – “SW-DOWN” Input from DOWN keypad. Grounding this pin will simulate the ‘down’ pad being pressed.

Pin 4 - Ground – common to one side of each of the keypads.

J4 – Input from Key Switch

For reference, wiring diagram 33530H37 provides a schematic depiction of the key switch and **J4**. Pin 1 is toward the top of the board, away from “J4” label.

Pin 1 – “+5V” Output to key switch

Pin 2 – no label Enables battery charging circuit. When key is in position 2 (/), +5 VDC will be present. +4.3ish VDC (diode drop from +5) will be present on this pin when the key switch is in position 3.

Pin 3 – “ALARM-EN” Signal to CPU board that the key switch is in position 3 (-). Also enables battery charging circuit via a diode to Pin 2. +5 VDC will be present on this pin when active.

Pin 4 – “PS15” Primary DC voltage supplied to key switch. Typically, this will measure 11.7ish VDC. During normal operation, as long as line voltage is supplied to the unit, this should always be present. See additional information below.

Pin 5 & 6 – “PS15R” This is the signal that initiates a power up. When the key is in position 2 or 3 (/ or -), the DC voltage measured at pin 1 should be present here. See additional information on page 6.

J5 - Serial Communications in from CPU board

The signals at J5 can not be accurately measured with a typical DMM. Some activity maybe discerned with a DMM but an oscilloscope is really required to capture these signals. Plus, measuring points on the 326301H01 are essentially inaccessible. Troubleshooting is basically limited to swapping harnesses or assemblies.

Pin 1 – Serial Data


Pin 2 – Latch

Pin 3 – Blank


Pin 4 – Ground

Pin 5 – Serial Clock

Pin 6 – Not used.



326276H01 Display Board Ladder Diagram




Cardinal Control Panel Overlays



LED Blue, 2010 design

314720H01 = +4/5° Glass door units with push button light switch.

314720H02 = +5° Solid door and all -20°/-30° units.



Black, 2018 design


327745H04 = +5° Solid door and all -20°/-30° units.

327745H05 = +4/5° Glass door units with integrated light switch.

TSX Thermo Scientific Overlays



327745H01 = +4/5° Glass door units with push button light switch.




327745H02 = +5° Solid door and all -20°/-30° units.

327745H03 = +4/5° Glass door units with integrated light switch.

H02 and H03 includes Dual Language English and French

Control Panel Information



The control panel is located on the top right side of your freezer. You can use the three pushbuttons (#5, #8, and #9 in *Figure above*) to change the temperature display (#1) or to adjust temperature and alarm setpoints as given in. The thermometer display (#2) provides a quick visual indicator of current cabinet temperature and alarm conditions.

1. **Main temperature display** - during normal operation, shows cabinet temperature in degrees Celsius, as measured by the primary sensor inside the cabinet. You can use the buttons to display other values such as setpoints and highest and lowest recorded temperatures. The number in the main display flashes when the value can be modified.
2. **Thermometer** - shows cabinet temperature and alarm conditions. There are 10 horizontal bars: 9 are displayed during normal operation, the tenth (top) bar indicates a warm alarm condition. The number of bars illuminated indicates approximate cabinet temperature. With the default settings, 4 to 6 bars illuminated indicate that the cabinet is at desired setpoint. For example, suppose that the cabinet temperature setpoint is -30°C and that the cold and warm alarm setpoints are -35°C and -20°C. Then the number of bars illuminated indicates cabinet temperature as follows:

| Refrigerator Bars Displayed | Temperature (°C) | Freezer Bars Displayed | Temperature (°C) |
|-----------------------------|------------------|------------------------|------------------|
| Bulb Only | 2 (Cold Alarm) | Bulb Only | -35 (Cold Alarm) |
| 1 Bar | 2.6 | 1 Bar | -33.5 |
| 2 Bars | 3.2 | 2 Bars | -32 |
| 3 Bars | 3.8 | 3 Bars | -30.5 |
| 4 Bars | 4.4 | 4 Bars | -29 |
| 5 Bars | 5 | 5 Bars | -27.5 |
| 6 Bars | 5.6 | 6 Bars | -26 |
| 7 Bars | 6.2 | 7 Bars | -24.5 |
| 8 Bars | 6.8 | 8 Bars | -23 |
| 9 Bars | 7.4 | 9 Bars | -21.5 |
| 10 Bars | 8 (Warm Alarm) | 10 Bars | -20 (Warm Alarm) |

When cabinet temperature exceeds the warm alarm setpoint, the top bar of the thermometer flashes. When temperature is lower than the cold alarm setpoint, the bulb flashes. When you are in programming mode the thermometer shows the setpoint value you are changing.

3. **Power failure** - illuminates when the main power supply is interrupted.

In this case the audible alarm also sounds.

4. **Service mode** - illuminates when the controller is in service programming mode.

5. **Increase** - pushbutton used to increase setpoint values in programming mode and for various display functions.

6. **Door ajar** - illuminates when the freezer door is open longer than the duration specified in the service parameters table; default is approximately 3 minutes (when the alarm is activated and the key switch is turned to the alarm position).

7. **Battery low** - illuminates when the backup battery is low

8. **Decrease** - pushbutton used to decrease setpoint values in programming mode and for various display functions.

9. **Scan** - pushbutton used to change the main display and for various other functions.

10. **Audible alarm** - illuminates during warm and cold alarm conditions.

11. **Key Switch** - switch used to turn the power and alarms on and off for the unit.

12. **Alarm On** - Setpoint Security - when the key switch is in the alarm on position, the alarms are on and the setpoints cannot be modified.

13. **Power On** - when the key switch is in the on position, the unit power is on with no alarms active.

14. **Off** - when the key switch is in the off position, the unit is off.

+4° / +5° Refrigerator Alarm Functions

| Alarm Description | Display Activity | Audible | Has Been | Remote Alarm Activated |
|---|--|--|--|------------------------|
| Warm Alarm | Flash uppermost bar of thermometer, flash speaker icon | * Yes * | Yes - slowly flash speaker icon and uppermost bar of thermometer alternately | Yes |
| Cold Alarm | Flash thermometer bulb, flash speaker icon | * Yes * | Yes - slowly flash speaker icon and thermometer bulb alternately | Yes |
| Power Failure | Reduced display intensity, lightning bolt on (no flash), flash temperature display approximately every 3 seconds | * Yes * | Yes - slowly flash speaker icon and lightning bolt alternately | Yes |
| Battery Low | Battery Icon On (no flash) | No | No | No |
| Door Ajar | Door Icon On (no flash) | No-if less than ~3 minutes, *Yes*-if more than 3 minutes | No | No |
| | | | | |
| | | | | |
| * Mute will silence audible alarm, with a ~10 minute ring back interval, display will remain flashing. If key switch is moved from position 3 to position 2, alarms will be disabled including remote alarms: audible and visuals will go away. | | | | |
| | | | | |
| | | | | |
| | | | | |
| "Err" occurs when the CPU cannot read the input channel - applies to sensors - see codes in System level troubleshooting for complete list | | | | |
| | | | | |
| | | | | |
| To clear Has Been Alarms, press Up and Down buttons simultaneously | | | | |

-20° / -30° Freezer Alarm Functions

| Alarm Description | Display Activity | Audible | Has Been | Remote Alarm Activated |
|---|--|--|---|------------------------|
| Warm Alarm | Flash upper most bar of thermometer, flash speaker icon | * Yes * | Yes - slowly flash speaker icon and upper most bar of thermometer alternately | Yes |
| Cold Alarm | Flash thermometer bulb, flash speaker icon | * Yes * | Yes - slowly flash speaker icon and thermometer bulb alternately | Yes |
| Power Failure | Reduced display intensity, lightning bolt on (no flash), flash temperature display approximately every 3 seconds | * Yes * | Yes - slowly flash speaker icon and lightning bolt alternately | Yes |
| Battery Low | Battery Icon On (no flash) | No | No | No |
| Door Ajar | Door Icon On (no flash) | No - if less than ~3 minutes. *Yes* - if more than 3 minutes | No | No |
| | | | | |
| * Mute will silence audible alarm, with a ~10 minute ring back interval, display will remain flashing. If key switch is moved from position 3 to position 2, alarms will be disabled including remote alarms: audible and visuals will go away. | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| "Err" occurs when the CPU cannot read the input channel - applies to sensors - see codes in System level troubleshooting for complete list | | | | |
| | | | | |
| | | | | |
| | To clear Has Been Alarms, press Up and Down buttons simultaneously | | | |

Feature & Alarm Descriptions

Refrigerator / Blood Bank / Freezer

| Feature | Complete set of Features/Alarms |
|------------------------------|---------------------------------|
| Temperature Display | xxx |
| Three Button Membrane | xxx |
| Cold/Warm Excursions | xxx |
| 3 position Key switch | xxx |
| Light Switch (if applicable) | xxx |
| Power Failure Alarm | xxx |
| High Temperature Alarm | xxx |
| Low Temperature Alarm | xxx |
| Door Ajar Alarm | xxx |
| Audible Alarm | xxx |
| Display Low Bottle (BB only) | xxx |
| Remote Alarm Contacts | xxx |
| Alarm Test | xxx |
| Battery Low Alarm | xxx |
| Thermometer Gauge | xxx |
| Setpoint Security | xxx |

Alarm Test

Alarm test is invoked by pressing the square button and increment at the same time:

- 4°, 5° units exercise the warm alarm first, then the cold alarm.
 - -20° & -30° units exercise the warm alarm.
 - If the alarm test is not successful within approximately 10 minutes after it is initiated, the service icon will be illuminated and the test is terminated.
 - The excursion feature is disabled from the time that an alarm test is initiated until approximately 5 minutes after it is terminated.
 - Unit returns to normal operation after termination whether test passes or fails.
-

Audible Alarm

Occurs for active warm alarm, cold alarm, and power failure. Pressing the square button will mute an audible alarm. If the alarm remains active, ring back occurs after approximately 10 minutes.

An audible alarm will occur if the door remains ajar for approximately 3 minutes, ring back occurs after approximately 10 minutes.

Battery Low Alarm

Battery voltage level is checked 4 hours after initial start up and every 8 hours thereafter. If the battery voltage level is found to be below approximately 9.8 volts DC, the battery icon will be illuminated.

NOTE: If the battery icon is illuminated and the battery is replaced, the icon will remain illuminated until the next test is performed (worst case 8 hours). Test for battery low cannot be performed by disconnecting the battery from the circuit, it must be tested with a discharged battery.

Estimated support time for the 3.3aH battery is 12 hours.

Cold Temperature Alarm

Alarm is invoked if the cabinet temperature (offset included), as measured in the upper bottle, exceeds the cold alarm setpoint. If this occurs, the following takes place:

- An audible alarm is invoked.
 - The bulb of the thermometer, the displayed cabinet temperature, and the speaker icon flash every second.
 - If the alarm condition disappears, the prior alarm will be displayed as a 'has been' alarm: the cabinet temperature will no longer flash, the bulb and the speaker icon will flash every other second and be out of phase from one another.
-

Display Low Bottle (Blood Banks only)

Lower bottle temperature is displayed by pressing the square button and decrement at the same time. When the buttons are released, the temperature displayed will default back to the upper bottle. Non Blood Banks show default "Err".

Door Ajar Alarm

Alarm is invoked if the door switch closes (door opens). If this occurs, the following takes place:

- The door icon is illuminated.
- The evaporator fans are disabled for 6 minutes only, then restart even if the door is not closed (except -20° units).
- An audible alarm will occur if the door remains open for approximately 3 minutes.

***NOTE – 45' sliders do not have door switches**

Power Failure Alarm

Feature is based on the detection of line frequency. If frequency is not detected, the following will take place:

- NVRAM is updated.
 - All relays are turned off.
 - Power Failure Icon is illuminated, if other alarm conditions are occurring, the associated icon will be illuminated. Bars of the thermometer are still active.
NOTE: If a 'has been' alarm was occurring when the power failure takes place, the power failure alarm will supercede the 'has been' alarm (this may result in the top thermometer bar or bulb remaining illuminated during the alarm).
 - An audible alarm is invoked.
 - The display intensity is reduced by ~75%.
 - The cabinet temperature is displayed every ~ three seconds.
NOTE: user interface will still operate in a power failure mode: excursions, setpoints, etc. are still accessible.
 - When line frequency is restored, the display will return to normal operation. The power failure icon and the speaker icon will flash every other second to convey a 'has been' condition.
 - **The compressor will not restart for ~ 7 minutes.**
-


Warm Temperature Alarm

Alarm is invoked if the cabinet temperature (offset included), as measured in the upper bottle, exceeds the warm alarm setpoint. If this occurs, the following takes place:

- An audible alarm is invoked.
 - The top bar of the thermometer, the displayed cabinet temperature, and the speaker icon flash every second.
 - If the alarm condition disappears, the prior alarm will be displayed as a 'has been' alarm: the cabinet temperature will no longer flash, the top bar and the speaker icon will flash every other second and be out of phase from one another.
-

Remote Alarm Contacts

Only enabled while the key switch is in the third position. Contacts will switch states during an active temperature alarm or power failure alarm.



To install the remote alarm, make the following connections:

1. Connect the COMMON terminal on the cabinet switch to the COMMON wire on the alarm.
- 2a. To get an alarm when the switch contacts open, connect the OPEN ON FAIL terminal on the cabinet to the OPEN ON FAIL wire on the alarm.
- 2b. To get an alarm when the switch contacts close, connect the CLOSE ON FAIL terminal on the cabinet to the CLOSE ON FAIL wire on the alarm. The COMMON and CLOSE ON FAIL wires must be tied together in this application.
3. Plug the alarm system service cord into an electrical outlet.
The contacts will trip in the event of a power outage, high temperature alarm or low temperature alarm.

- All units have factory-installed local alarm contacts that can be used for remote alarm systems.
- The maximum distance between a refrigerator and a remote alarm depends on the wire gauge used.
- Remote alarm terminals are located at the rear of the machine compartment. The three terminals are: COMMON, OPEN ON FAIL (Normally Closed), and CLOSE ON FAIL (Normally Open)

| Wire Gauge | Total Wire Length (feet) | Distance to Alarm 1/2 Wire Length (feet) |
|------------|--------------------------|--|
| 20 | 530 | 265 |
| 18 | 840 | 420 |
| 16 | 1330 | 665 |
| 14 | 2120 | 1060 |
| 12 | 3370 | 1685 |

Setpoint Security

When the key is turned to the alarm position, the parameters are visible, but not adjustable. Parameters must be flashing to be settable (Alarm Off mode).

Thermometer Gauge

The number of bars illuminated within the thermometer is indicative of the cabinet temperature relative to the control setpoints. 4°/5° units and -20°/-30° are fundamentally different in this area of operation. The number of bars illuminated is determined in the following manner:

- 4° units

$$\# \text{ bars illuminated} = \frac{[\text{Cabinet Temperature} - \text{Low Alarm Setpoint}]}{[\text{Warm Alarm Setpoint} - \text{Low Alarm Setpoint}]} * 10 \text{ (bars)}$$

- -20° & -30° units








If Cabinet Temperature < Setpoint,

$$\# \text{ bars illuminated} = \frac{[\text{Cabinet Temperature} - \text{Low Alarm Setpoint}]}{[\text{Setpoint} - \text{Low Alarm Setpoint}]} * 5 \text{ (bars)}$$

If Cabinet Temperature > Setpoint,

$$\# \text{ bars illuminated} = \frac{[\text{Cabinet Temperature} - \text{Setpoint}]}{[\text{Warm Alarm Setpoint} - \text{Setpoint}]} * 5 \text{ (bars)} + 5 \text{ (bars)}$$




Basic Functions

| Function | Meaning | Sequence | Display |
|----------------------------|--|---|--|
| Normal operation | Default display while freezer is running. | | Temperature display and control panel thermometer icon show cabinet temperature. |
| Coldest logged temperature | Show coldest cabinet temperature since last startup or reset. | Press  | Display shows coldest logged temperature while button is pressed. |
| Warmest logged temperature | Show warmest cabinet temperature since last startup or reset. | Press  | Display shows warmest logged temperature while button is pressed. |
| Mute | Silence audible alarm for approximately 6 minutes. | Press  | Display and thermometer show cabinet temperature, alarm icon continues to flash. |
| Reset | Return to default display after alarm condition, clears temperature log. | Press and hold  and  simultaneously. | Excursion values are reset; temperature display shows cabinet temperature. |
| Alarm Test | Test by simulating warm alarm. Key switch must be in alarm mode. | Press  and  simultaneously, hold for approximately 5 seconds. | First, the display will show “AtSt” to show entry of the test. Display and thermometer show simulated cabinet temperatures (warm), alarms flash and sound as appropriate. Alarms clear when test is completed. |

Programming Mode

You can enter the programming mode by pressing the Scan button and holding for approximately 5 seconds. The display will then flash “Prg” to indicate that you have entered the programming mode. Use (UP) and (Down) buttons to modify the values and press the scan button to save the value and go to the next screen.

Note: If the alarms are enabled, all parameters are read only and cannot be edited. Programming Mode will time out after 30 seconds of inactivity.

| Screen/Button | | Function | Summary |
|---------------|--|------------------------|--|
| 1 | Initial Screen | Unit Setpoint | The temperature that the unit is set to run |
| 2 | Press Scan () | Cold Alarm Temperature | The temperature that will actuate the cold alarm (if active) |
| 3 | Press Scan () | Warm Alarm Temperature | The temperature that will actuate the warm alarm (if active) |
| 4 | Press Scan () | Exit Program Mode | Returns to normal operating screen |

Service Parameters






You can enter the service mode by pressing the Scan button and holding for approximately 10 seconds while in Programming mode. The display will then flash “SEr” to indicate that you have entered the service mode followed by the software checksum values flashing on the screen. The service icon will also illuminate. Use up, and down buttons to modify the values and press the scan button to save the value and go to the next screen.
















Note: If the alarms are enabled, the unit will not enter Service Mode.










Resetting any of the following parameter values could adversely affect the performance of your unit. Be sure to understand your product requirements before making any adjustment to the service parameter values. These settings very rarely need to be changed from normal operation. Call technical support if you have any questions prior to making any adjustments to service parameter values.









+4° / +5° Service Menu Functions

| Screen/Button | | Display | Function | Summary |
|---------------|--|----------|-------------------------------|--|
| 1 | Initial Screen | SEr 05 | Temp of unit | Allows the user to see what temperature the unit is designed to run. |
| | | Varies | CPU Checksum | Shows the current checksum of the CPU board for validation. |
| | | Varies | Relay Board Checksum | Shows the current checksum of the Relay board for validation. |
| | | tyP | Model Type | Unit type 05 for +5 refrigerator. This parameter should never be changed. |
| 2 | Press Scan () | CFg | Type of unit | This parameter should only be changed by trained service personnel. |
| 3 | Press Scan () | CuFt | Size of the unit | Setting for size of the unit. |
| 4 | Press Scan () | Pd oFSt | Control Probe Pulldown Offset | The difference in temperature between the control probe value and the average compartment temperature that determines when the control system changes from pull down mode to steady state operation. This can be useful to adjust if the average compartment temperature has an undershoot or overshoot relative to setpoint during pull down mode. This offset is also used on some refrigerator models as a response to door openings. |
| 5 | Press Scan () | Cnt oFSt | Control Probe Offset | The difference in temperature between the control probe value and the average compartment temperature that is used to control the cooling system to the unit setpoint during steady state operation. This can be useful to adjust if there is a discrepancy between the unit setpoint and the average cabinet temperature. |
| 6 | Press Scan () | dIS oFSt | Display Probe Offset | The difference in temperature between the display probe value seen on the User Interface display and the average compartment temperature at the unit setpoint during steady state operation. This can be useful to adjust the display during unit calibration procedures. |








| Screen/Button | | Display | Function | Summary |
|---------------|--|-----------|-----------------------------------|---|
| 7 | Press Scan () | Ctr CAI | Display Control Probe with Offset | (If present) Normally OFF (00). If ON (01), the Upper Bottle Temperature output on main display is replaced by the control probe temperature value with offset when exiting the service menu. |
| 8 | Press Scan () | Cnt ucl | Control Upper Hysteresis | The degree increase in temperature from the setpoint that will trigger the unit to begin cooling to prevent the unit from getting too warm. |
| 9 | Press Scan () | Cnt lcl | Control Lower Hysteresis | The degree decrease in temperature from the setpoint that will trigger the unit to stop cooling to prevent the unit from getting too cold. |
| 10 | Press Scan () | door AJAr | Door Alarm Time | The approximate time the door can remain open before the door ajar alarm activates in minutes. |
| 11 | Press Scan () | dEF Int | Defrost Interval | The time interval between defrost cycles in hours. |
| 12 | Press Scan () | dEF dur | Defrost Duration | The maximum amount of time the defrost cycle will run in minutes. |
| 13 | Press Scan () | dEF Hi | Defrost Temp Out | The maximum temperature that the defrost probe can reach before stopping the defrost cycle. |
| 14 | Press Scan () | qUA | Quality Mode | Normally OFF(00). Used for production line testing. |
| 15 | Press Scan () | bot | BOT Mode | Normally OFF(00). If ON(01), the unit is forced to run at coldest possible temperature for approximately 24 hours. Defrosts will still occur as scheduled. |
| 16 | Press Scan () | PEr | Perimeter Heater Duty Cycle | Controls the perimeter heater ON time (only applicable on swinging door models). |
| 17 | Press Scan () | Cnt Prb | Control Probe Temperature | Displays control probe temperature (Read only). |
| 18 | Press Scan () | Anb Prb | Ambient Probe Temperature | Displays ambient probe temperature in deck (Read only). |
| 19 | Press Scan () | d lo Prb | Lower Display Probe | Displays Err since this refrigerator does not have a lower display probe. |
| 20 | Press Scan () | dEF Prb | Defrost Probe Temperature | Displays defrost probe temperature (Read only). |
| 21 | Press Scan () | --- | Enter Defrost Cycle | Pressing the Scan button after displaying the defrost probe temperature while not in an alarm mode will force the unit into a defrost cycle. To avoid defrost cycle, leave unit screen idle for approximately 30 seconds. |















-20° Service Menu Functions

| Screen/Button | | Display | Function | Summary |
|---------------|--|----------|-----------------------------------|---|
| 1 | Initial Screen | SER -20 | Temp of unit | Allows the user to see what temperature the unit is designed to run. |
| | | Varies | CPU Checksum | Shows the current checksum of the CPU board for validation. |
| | | Varies | Relay Board Checksum | Shows the current checksum of the Relay board for validation. |
| | | tyP | Model Type | Unit type 20 for -20 Freezer. This parameter should never be changed. |
| 2 | Press Scan  | CFg | Type of unit | This parameter should only be changed by trained service personnel. |
| 3 | Press Scan  | CuFt | Size of the unit | Setting for Size of the unit. |
| 4 | Press Scan  | Pd oFSt | Control Probe Pulldown Offset | The difference in temperature between the control probe value and the average compartment temperature that determines when the control system changes from pull down mode to steady state operation. This can be useful to adjust if the average compartment temperature has an undershoot or overshoot relative to setpoint during pull down mode. This offset is also used on some freezer models as a response to door openings. |
| 5 | Press Scan  | Cnt oFSt | Control Probe Offset | The difference in temperature between the control probe value and the average compartment temperature that is used to control the cooling system to the unit setpoint during steady state operation. This can be useful to adjust if there is a discrepancy between the unit setpoint and the average cabinet temperature. |
| 6 | Press Scan  | dIS oFSt | Display Probe Offset | The difference in temperature between the display probe value seen on the User Interface display and the average compartment temperature at the unit setpoint during steady state operation. This can be useful to adjust the display during unit calibration procedures. |
| 7 | Press Scan  | Ctr CAI | Display Control Probe with Offset | (If present) Normally OFF (00). If ON (01), the Upper Bottle Temperature output on main display is replaced by the control probe temperature value with offset when exiting the service menu. |
| 8 | Press Scan  | Cnt ucl | Control Upper Hysteresis | The degree increase in temperature from the setpoint that will trigger the unit to begin cooling to prevent the unit from getting too warm. |

| Screen/Button | | Display | Function | Summary |
|---------------|--|-----------|-----------------------------|--|
| 9 | Press Scan () | Cnt lcl | Control Lower Hysteresis | The degree decrease in temperature from the setpoint that will trigger the unit to stop cooling to prevent the unit from getting too cold. |
| 10 | Press Scan () | door AJAr | Door Alarm Time | The approximate time the door can remain open before the door ajar alarm activates in minutes. |
| 11 | Press Scan () | qUA | Quality Mode | Normally OFF(00). Used for production line testing. |
| 12 | Press Scan () | bot | BOT Mode | Normally OFF(00). If ON(01), the unit is forced to run at coldest possible temperature for approximately 24 hours. |
| 13 | Press Scan () | PEr | Perimeter Heater Duty Cycle | Controls the perimeter heater ON time. |
| 14 | Press Scan () | Cnt Prb | Control Probe Temperature | Displays control probe temperature (Read only). |
| 15 | Press Scan () | Anb Prb | Ambient Probe Temperature | Displays ambient probe temperature in deck (Read only). |
| 16 | Press Scan () | --- | Exit Service Mode | Returns to normal operating screen. |

-30° Service Menu Functions

| Screen/Button | | Display | Function | Summary |
|---------------|--|----------|-----------------------------------|---|
| 1 | Initial Screen | SEr -30 | Temp of unit | Allows the user to see what temperature the unit is designed to run. |
| | | Varies | CPU Checksum | Shows the current checksum of the CPU board for validation. |
| | | Varies | Relay Board Checksum | Shows the current checksum of the Relay board for validation. |
| | | tyP | Model Type | Unit type 30 for -30 Freezer. This parameter should never be changed. |
| 2 | Press Scan () | CFg | Type of unit | This parameter should only be changed by trained service personnel. |
| 3 | Press Scan () | CuFt | Size of the unit | Setting for Size of the unit. |
| 4 | Press Scan () | Pd oFSt | Control Probe Pulldown Offset | The difference in temperature between the control probe value and the average compartment temperature that determines when the control system changes from pull down mode to steady state operation. This can be useful to adjust if the average compartment temperature has an undershoot or overshoot relative to setpoint during pull down mode. This offset is also used on some freezer models as a response to door openings. |
| 5 | Press Scan () | Cnt oFSt | Control Probe Offset | The difference in temperature between the control probe value and the average compartment temperature that is used to control the cooling system to the unit setpoint during steady state operation. This can be useful to adjust if there is a discrepancy between the unit setpoint and the average cabinet temperature. |
| 6 | Press Scan () | dIS oFSt | Display Probe Offset | The difference in temperature between the display probe value seen on the User Interface display and the average compartment temperature at the unit setpoint during steady state operation. This can be useful to adjust the display during unit calibration procedures. |
| 7 | Press Scan () | Ctr CAI | Display Control Probe with Offset | (If present) Normally OFF (00). If ON (01), the Upper Bottle Temperature output on main display is replaced by the control probe temperature value with offset when exiting the service menu. |
| 8 | Press Scan () | Cnt ucl | Control Upper Hysteresis | The degree increase in temperature from the setpoint that will trigger the unit to begin cooling to prevent the unit from getting too warm. |

| Screen/Button | | Display | Function | Summary |
|---------------|--|---------|------------------------------------|---|
| 8 | Press Scan  | Cnt lcl | Control Lower Hysteresis | The degree decrease in temperature from the set point that will trigger the unit to stop cooling to prevent the unit from getting too cold. |
| 9 | Press Scan  | AdP dEF | Adaptive Defrost | Setting for adaptive defrost ON (01) or OFF (00). |
| 10 | Press Scan  | dEF Int | Defrost Interval | The time interval between defrost cycles in hours. |
| 11 | Press Scan  | dEF dur | Defrost Duration | The maximum amount of the time the defrost cycle will run in minutes. Value will be "nA" for units with hot gas defrost |
| 12 | Press Scan  | dEF Hi | Defrost Temp Out | The maximum temperature that the defrost probe can reach before stopping the defrost cycle. |
| 13 | Press Scan  | PAn Hi | Drip Pan Temp Out | Maximum temperature drip pan probe can reach before stopping the defrost drip time. |
| 14 | Press Scan  | qUA | Quality Mode | Normally OFF(00). Used for production line testing. |
| 15 | Press Scan  | bot | BOT mode | Normally OFF(00). If ON(01), the unit is forced to run at coldest possible temperature for approximately 24 hours. Defrosts will still occur at scheduled. |
| 16 | Press Scan  | PEr | Perimeter heater Duty Cycle | Controls the perimeter heater ON time. |
| 17 | Press Scan  | Cnt prb | Control Probe Temperature | Displays control probe temperature (Read Only). |
| 18 | Press Scan  | Anb prb | Ambient Probe Temperature | Displays ambient probe temperature in deck (Read Only). |
| 19 | Press Scan  | dEF Prb | Defrost Probe Temperature | Displays defrost probe temperature (Read Only). |
| 20 | Press Scan  | PAn Prb | Defrost Drip Pan Probe Temperature | Displays defrost drip pan probe temperature. |
| 21 | Press Scan  | --- | Enter Defrost Cycle | Pressing the scan button after displaying the defrost Drip Pan probe temperature while not in an alarm mode will force the unit into a defrost cycle. To avoid defrost cycle, leave idle screen for approximately 30 seconds. |

Firmware History

| Vintage | EPROM version | EPROM Checksum | Relay Checksum | Approximate Production Dates | Notes |
|-----------------|---------------|----------------|----------------|--|---|
| Suffix 16 | G52B | 1AB1 | 4CE4 | Nov 2016 to May 2018 | 2330 Electric Defrost |
| Suffix 17 | G54A | 520E | 3CA0 | May 2018 to Apr 2019 | Switch to Hot Gas Defrost, H02 Board Not Working on Rev A units. |
| Suffix 17 | G54B | CF1B | DB6A | Apr 2019 to May 2019 | Allow new H02 board to work on REVA units |
| Suffix 17/18 | G54C | A678 | 9AC9 | May 2019 to Feb 2020 | Product Launch, Issue with ADP = 00 No Defrost on -30 |
| Suffix 17/18/19 | G54D | A901 | D284 | Feb 2020 to Nov 2020 | Fixed ADP=00 Issue , added Door Ajar, Moved control sensor on CW's, so suffix 19 - Config 03 |
| Suffix 17/19 | G54D | A901 | 3429 | Dec 2020 to May 2021 | Production fix for 2330 only |
| Suffix 17-20 | G54E | 36DE | 9862 | May 2021 to current | Indroduction of the 2320HA FMS model; 2320 standerize MAIA inverter 3020 Cond. unit, no Transformer; 2330 new compressor/inverter. 2304-5 Standardization |
| Icing Fix | | | | | |
| Vintage | EPROM version | EPROM Checksum | Relay Checksum | Approximate Production Dates | Notes |
| Suffix 17 | G57A | DC38 | F0DD | TSX 5030 Icing fix | Dec 2019, fix for 5030 units prior to 11-8-2019 Fix for Adaptive defrost of 00 equals no defrost at all ; also added Door Ajar |
| Suffix 17 | G57B | 1094 | FE0C | Nov. 2020 fix for 1230, 2330, and 3030 | Nov 2020, fix for 1230, 2330, and 3030 in addition to 5030. Added Control Calibration (Ctr Cal) to the service menu |

The new user manual states that the Type of Unit for the Config Parameter should not be changed. That is correct when the unit is produced but may not apply for service boards.

Cfg = 00 is used for Suffix 16 models that were produced from 2016 to May 2018.

Cfg = 01 is for most 2018 and 2019 models.

Cfg = 02 is for the 2019 version of the coldwall units that had the addition of a display bottle and a PEP.

Cfg = 03 is for the 2020 version of the coldwall units that had a longer red control probe mounted to the lower left side wall with ULT sensor bracket.

Cfg = 04 for all units May 2021 to current; changes for the 2320 units with 3020 style condensing unit, and the addition of the new service parameter Control Calibration.

The 2019 menus now include offsets that were previously hidden and hard coded. Pull down offset and Control probe offset have default values based on cabinet size and temperature, but are now user settable for customers who have their own validation protocols that needed access to those settings.

Parameter Defaults by Eprom Revision and Suffix

| Suffix 17 G54 Rev. A Paramters May 2018 - Apr. 2019 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
|---|----------------------|-----------|----------------|------------|----------------------|-------------------|------------------|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|----------------------|-----------------|------------------|------------------|---------------------|-----------------|-----------------|----------------------|-----------------|-----------------|---------------------|-----------|----------|---------|---------|
| "+4 Refrigerator" | | | | | | "+5 Refrigerator" | | | | | | | | | | | | | | | | | | | | | | | |
| 4 | System Parameter | Units | Display "text" | Value | Range | 5 | System Parameter | Units | Display "text" | Value | Range | | | | | | | | | | | | | | | | | | |
| PROG MODE | Set point | °C | na | 4.0 | 3 to 7 | PROG MODE | Set point | °C | na | 5.0 | 3 to 7 | | | | | | | | | | | | | | | | | | |
| | Warm Alarm | °C | na | 5.5 | (sp+1) to 15 | | PROG MODE | Warm Alarm | °C | na | 8.0 | (sp+1) to 15 | | | | | | | | | | | | | | | | | |
| | Cold Alarm | °C | na | 2.0 | 0 to (sp-1) | | | PROG MODE | Cold Alarm | °C | na | 2.0 | 0 to (sp-1) | | | | | | | | | | | | | | | | |
| SERVICE MODE | "SER" | na | SEr | na | na | SERVICE MODE | | | "SER" | na | SEr | na | na | | | | | | | | | | | | | | | | |
| | Temp of Unit | °C | na | 4 | 4,20,30,5 | | SERVICE MODE | | Temp of Unit | °C | na | 5 | 4,20,30,5 | | | | | | | | | | | | | | | | |
| | CPU Checksum | base 16 # | na | 520E | na | | | SERVICE MODE | CPU Checksum | base 16 # | na | 520E | na | | | | | | | | | | | | | | | | |
| | Relay Checksum | base 16 # | na | 3CA0 | na | | | | SERVICE MODE | Relay Checksum | base 16 # | na | 3CA0 | na | | | | | | | | | | | | | | | |
| | Model Type | coded # | tyP | 4 | 4, 20, 30, 5, 55, 66 | | | | | SERVICE MODE | Model Type | coded # | tyP | 5 | 4, 20, 30, 5, 55, 66 | | | | | | | | | | | | | | |
| | Configuration | coded # | CFg | 01 | 00, 01 | | | | | | SERVICE MODE | Configuration | coded # | CFg | 01 | 00, 01 | | | | | | | | | | | | | |
| | Cubic Foot | cu ft | CuFt | (cab size) | 23, 50 | | | | | | | SERVICE MODE | Cubic Foot | cu ft | CuFt | (cab size) | 23, 45, 50 | | | | | | | | | | | | |
| | Display Offset | °C | oFSt | 0.0 | +/- 3.0 | | | | | | | | SERVICE MODE | Display Offset | °C | oFSt | 0.0 | +/- 3.0 | | | | | | | | | | | |
| | Control UCL | °C | Cnt ucl | 0.8 | 0 to 5 | | | | | | | | | SERVICE MODE | Control UCL | °C | Cnt ucl | 0.8 | 0 to 5 | | | | | | | | | | |
| | Control LCL | °C | Cnt lcl | -0.2 | -5 to 0 | | | | | | | | | | SERVICE MODE | Control LCL | °C | Cnt lcl | -0.2 | -5 to 0 | | | | | | | | | |
| | Defrost Interval | hours | dEF Int | 1 | 1 to 12 | | | | | | | | | | | SERVICE MODE | Defrost Interval | hours | dEF Int | 1 | 1 to 12 | | | | | | | | |
| | Defrost Duration | minutes | dEF dur | 15 | 15 to 30 | | | | | | | | | | | | SERVICE MODE | Defrost Duration | minutes | dEF dur | 15 | 15 to 30 | | | | | | | |
| | Defr. Hi (temp out) | °C | dEF Hi | 5 | 1 to 15 | | | | | | | | | | | | | SERVICE MODE | Defr. Hi (temp out) | °C | dEF Hi | 5 | 1 to 15 | | | | | | |
| | Quality Mode | on/off | qUA | 00 | 00, 01 | | | | | | | | | | | | | | SERVICE MODE | Quality Mode | on/off | qUA | 00 | 00, 01 | | | | | |
| | BOT Request | on/off | bot | 00 | 00, 01 | | | | | | | | | | | | | | | SERVICE MODE | BOT Request | on/off | bot | 00 | 00, 01 | | | | |
| | Perimeter Heater PWM | % | Per | 35 | 20 to 100 | | | | | | | | | | | | | | | | SERVICE MODE | Perimeter Heater PWM | % | Per | 35 | 20 to 100 | | | |
| | Cntr.Probe | °C | Cnt Prb | °C | reading | | | | | | | | | | | | | | | | | SERVICE MODE | Cntr.Probe | °C | Cnt Prb | °C | reading | | |
| | Amb. Probe | °C | Anb Prb | °C | reading | | | | | | | | | | | | | | | | | | SERVICE MODE | Amb. Probe | °C | Anb Prb | °C | reading | |
| | Lower Display Probe | °C | d lo Prb | °C | reading | | | | | | | | | | | | | | | | | | | SERVICE MODE | Lower Display Probe | °C | d lo Prb | °C | reading |
| | Defr.Prb | °C | dEF Prb | °C | reading | | | | | | | | | | | | | | | | | | | | SERVICE MODE | Defr.Prb | °C | dEF Prb | °C |

| Suffix 17 G54 Rev. A Parameters May 2018 - Apr. 2019 | | | | | | | | | | | | |
|--|----------------------|-----------|----------------|-------|----------------------|---------------------|----------------------|------------------|----------------|----------|---------------|----------------------|
| "-20 Freezer" | | | | | | "-30 Freezer" | | | | | | |
| 20 | System Parameter | Units | Display "text" | Value | Range | 30 | System Parameter | Units | Display "text" | Value | Range | |
| PROG MODE | Set point | °C | na | -20.0 | -25 to -15 | PROG MODE | Set point | °C | na | -30.0 | -35 to -15 | |
| | Warm Alarm | °C | na | -10.0 | (sp+5) to -10 | | Warm Alarm | °C | na | -20.0 | (sp+5) to -10 | |
| | Cold Alarm | °C | na | -25.0 | -35 to (sp-5) | | | Cold Alarm | °C | na | -35.0 | -40 to (sp-5) |
| SERVICE MODE | "SER" | na | SER | na | na | SERVICE MODE | | | "SER" | na | SER | na |
| | Temp of Unit | °C | na | 20 | 4,20,30,5 | | Temp of Unit | | °C | na | 30 | 4,20,30,5 |
| | CPU Checksum | base 16 # | na | 520E | na | | | CPU Checksum | base 16 # | na | 520E | na |
| | Relay Checksum | base 16 # | na | 3CA0 | na | Relay Checksum | | | base 16 # | na | 3CA0 | na |
| | Model Type | coded # | tyP | 20 | 4, 20, 30, 5, 55, 66 | | Model Type | | coded # | tyP | 30 | 4, 20, 30, 5, 55, 66 |
| | Configuration | coded # | CFg | 01 | 00, 01 | | | Configuration | coded # | CFg | 01 | 00, 01 |
| | Cubic Foot | cu ft | CuFt | 23 | 23 | Cubic Foot | | | cu ft | CuFt | (cab size) | 23, 50 |
| | Display Offset | °C | oFSt | 0.0 | +/- 10 | | Display Offset | | °C | oFSt | 0.0 | +/- 10 |
| | Control UCL | °C | Cnt ucl | 0.5 | 0 to 5 | | | Control UCL | °C | Cnt ucl | 0.5 | 0 to 5 |
| | Control LCL | °C | Cnt lcl | -0.5 | -5 to 0 | Control LCL | | | °C | Cnt lcl | -0.5 | -5 to 0 |
| | Quality Mode | on/off | qUA | 00 | 00, 01 | | Defrost Interval | | hours | dEF Int | 6 | 4 to 8 |
| | BOT Request | on/off | bot | 00 | 00, 01 | | | Defrost Duration | minutes | dEF dur | 20 | 20 to 30 |
| | Perimeter Heater PWM | % | Per | 64 | 20 to 100 | Defr. Hi (temp out) | | | °C | dEF Hi | 1 | 1 to 10 |
| | Cntr.Probe | °C | Cnt Prb | °C | reading | | Drip Pan temp out | | °C | Pan Hi | 3 | 0 to 10 |
| | Amb. Probe | °C | Anb Prb | °C | reading | | | Quality Mode | on/off | qUA | 00 | 00, 01 |
| | | | | | | BOT Request | | | on/off | bot | 00 | 00, 01 |
| | | | | | | | Perimeter Heater PWM | | % | Per | 90 | 20 to 100 |
| | | | | | | | | Cntr.Probe | °C | Cnt Prb | °C | reading |
| | | | | | | Amb. Probe | | | °C | Anb Prb | °C | reading |
| | | | | | | | Lower Display Probe | | °C | d lo Prb | °C | reading |
| | | | | | | | | Defr.Prb | °C | dEF Prb | °C | reading |
| | | | | | | Drip Pan Probe | | | °C | PAn Prb | °C | reading |

Suffix 17 G54 Rev. B Paramters May 2018 - Apr. 2019

| "+4 Refrigerator" | | | | | | "+5 Refrigerator" | | | | | |
|-------------------|----------------------|-----------|----------------|------------|----------------------|-------------------|----------------------|-----------|----------------|------------|----------------------|
| 4 | System Parameter | Units | Display "text" | Value | Range | 5 | System Parameter | Units | Display "text" | Value | Range |
| PROG | Set point | °C | na | 4.0 | 3 to 7 | PROG | Set point | °C | na | 5.0 | 3 to 7 |
| MODE | Warm Alarm | °C | na | 5.5 | (sp+1) to 15 | MODE | Warm Alarm | °C | na | 8.0 | (sp+1) to 15 |
| | Cold Alarm | °C | na | 2.0 | 0 to (sp-1) | | Cold Alarm | °C | na | 2.0 | 0 to (sp-1) |
| | "SER" | na | SER | na | na | | "SER" | na | SER | na | na |
| SERVICE | Temp of Unit | °C | na | 4 | 4,20,30,5 | SERVICE | Temp of Unit | °C | na | 5 | 4,20,30,5 |
| MODE | CPU Checksum | base 16 # | na | CF1B | na | MODE | CPU Checksum | base 16 # | na | CF1B | na |
| | Relay Checksum | base 16 # | na | DB6A | na | | Relay Checksum | base 16 # | na | DB6A | na |
| | Model Type | coded # | tyP | 4 | 4, 20, 30, 5, 55, 66 | | Model Type | coded # | tyP | 5 | 4, 20, 30, 5, 55, 66 |
| | Configuration | coded # | CFg | 01 | 00, 01 | | Configuration | coded # | CFg | 01 | 00, 01 |
| | Cubic Foot | cu ft | CuFt | (cab size) | 23, 50 | | Cubic Foot | cu ft | CuFt | (cab size) | 23, 45, 50 |
| | Display Offset | °C | oFSt | 0.0 | +/- 3.0 | | Display Offset | °C | oFSt | 0.0 | +/- 3.0 |
| | Control UCL | °C | Cnt ucl | 0.8 | 0 to 5 | | Control UCL | °C | Cnt ucl | 0.8 | 0 to 5 |
| | Control LCL | °C | Cnt lcl | -0.2 | -5 to 0 | | Control LCL | °C | Cnt lcl | -0.2 | -5 to 0 |
| | Defrost Interval | hours | dEF Int | 1 | 1 to 12 | | Defrost Interval | hours | dEF Int | 1 | 1 to 12 |
| | Defrost Duration | minutes | dEF dur | 15 | 15 to 30 | | Defrost Duration | minutes | dEF dur | 15 | 15 to 30 |
| | Defr. Hi (temp out) | °C | dEF Hi | 5 | 1 to 15 | | Defr. Hi (temp out) | °C | dEF Hi | 5 | 1 to 15 |
| | Quality Mode | on/off | qUA | 00 | 00, 01 | | Quality Mode | on/off | qUA | 00 | 00, 01 |
| | BOT Request | on/off | bot | 00 | 00, 01 | | BOT Request | on/off | bot | 00 | 00, 01 |
| | Perimeter Heater PWM | % | Per | 35 | 20 to 100 | | Perimeter Heater PWM | % | Per | 35 | 20 to 100 |
| | Cntr.Probe | °C | Cnt Prb | °C | reading | | Cntr.Probe | °C | Cnt Prb | °C | reading |
| | Amb. Probe | °C | Anb Prb | °C | reading | | Amb. Probe | °C | Anb Prb | °C | reading |
| | Lower Display Probe | °C | d lo Prb | °C | reading | | Lower Display Probe | °C | d lo Prb | °C | reading |
| | Defr.Prb | °C | dEF Prb | °C | reading | | Defr.Prb | °C | dEF Prb | °C | reading |

| Suffix 17 G54 Rev. B Paramters May 2018 - Apr. 2019 | | | | | | | | | | | |
|---|----------------------|-----------|----------------|-------|----------------------|---------------|----------------------|-----------|----------------|------------|----------------------|
| "-20 Freezer" | | | | | | "-30 Freezer" | | | | | |
| 20 | System Parameter | Units | Display "text" | Value | Range | 30 | System Parameter | Units | Display "text" | Value | Range |
| PROG | Set point | °C | na | -20.0 | -25 to -15 | PROG | Set point | °C | na | -30.0 | -35 to -15 |
| MODE | Warm Alarm | °C | na | -10.0 | (sp+5) to -10 | MODE | Warm Alarm | °C | na | -20.0 | (sp+5) to -10 |
| | Cold Alarm | °C | na | -25.0 | -35 to (sp-5) | | Cold Alarm | °C | na | -35.0 | -40 to (sp-5) |
| | "SER" | na | SEr | na | na | | "SER" | na | SEr | na | na |
| SERVICE | Temp of Unit | °C | na | 20 | 4,20,30,5 | SERVICE | Temp of Unit | °C | na | 30 | 4,20,30,5 |
| MODE | CPU Checksum | base 16 # | na | CF1B | na | MODE | CPU Checksum | base 16 # | na | CF1B | na |
| | Relay Checksum | base 16 # | na | DB6A | na | | Relay Checksum | base 16 # | na | DB6A | na |
| | Model Type | coded # | tyP | 20 | 4, 20, 30, 5, 55, 66 | | Model Type | coded # | tyP | 30 | 4, 20, 30, 5, 55, 66 |
| | Configuration | coded # | CFg | 01 | 00, 01 | | Configuration | coded # | CFg | 01 | 00, 01 |
| | Cubic Foot | cu ft | CuFt | 23 | 23 | | Cubic Foot | cu ft | CuFt | (cab size) | 23, 50 |
| | Display Offset | °C | oFSt | 0.0 | +/- 10 | | Display Offset | °C | oFSt | 0.0 | +/- 10 |
| | Control UCL | °C | Cnt ucl | 0.5 | 0 to 5 | | Control UCL | °C | Cnt ucl | 0.5 | 0 to 5 |
| | Control LCL | °C | Cnt lcl | -0.5 | -5 to 0 | | Control LCL | °C | Cnt lcl | -0.5 | -5 to 0 |
| | Quality Mode | on/off | qUA | 00 | 00, 01 | | Defrost Interval | hours | dEF Int | 6 | 4 to 8 |
| | BOT Request | on/off | bot | 00 | 00, 01 | | Defrost Duration | minutes | dEF dur | 20 | 20 to 30 |
| | Perimeter Heater PWM | % | Per | 64 | 20 to 100 | | Defr. Hi (temp out) | °C | dEF Hi | 1 | 1 to 10 |
| | Cntr.Probe | °C | Cnt Prb | °C | reading | | Drip Pan temp out | °C | Pan Hi | 3 | 0 to 10 |
| | Amb. Probe | °C | Anb Prb | °C | reading | | Quality Mode | on/off | qUA | 00 | 00, 01 |
| | | | | | | | BOT Request | on/off | bot | 00 | 00, 01 |
| | | | | | | | Perimeter Heater PWM | % | Per | 90 | 20 to 100 |
| | | | | | | | Cntr.Probe | °C | Cnt Prb | °C | reading |
| | | | | | | | Amb. Probe | °C | Anb Prb | °C | reading |
| | | | | | | | Lower Display Probe | °C | d lo Prb | °C | reading |
| | | | | | | | Defr.Prb | °C | dEF Prb | °C | reading |
| | | | | | | | Drip Pan Probe | °C | PAn Prb | °C | reading |

Suffix 17 G54 Rev. C Paramters Apr. 2019 - Feb. 2020

| "+4 Refrigerator" | | | | | | "+5 Refrigerator" | | | | | |
|-------------------|----------------------|-----------|----------------|------------|----------------------|-------------------|----------------------|-----------|----------------|------------|----------------------|
| 4 | System Parameter | Units | Display "text" | Value | Range | 5 | System Parameter | Units | Display "text" | Value | Range |
| PROG | Set point | °C | na | 4.0 | 3 to 7 | PROG | Set point | °C | na | 5.0 | 3 to 7 |
| MODE | Warm Alarm | °C | na | 5.5 | (sp+1) to 15 | MODE | Warm Alarm | °C | na | 8.0 | (sp+1) to 15 |
| | Cold Alarm | °C | na | 2.0 | 0 to (sp-1) | | Cold Alarm | °C | na | 2.0 | 0 to (sp-1) |
| | "SER" | na | SEr | na | na | | "SER" | na | SEr | na | na |
| SERVICE | Temp of Unit | °C | na | 4 | 4,20,30,5 | SERVICE | Temp of Unit | °C | na | 5 | 4,20,30,5 |
| MODE | CPU Checksum | base 16 # | na | A678 | na | MODE | CPU Checksum | base 16 # | na | A678 | na |
| | Relay Checksum | base 16 # | na | 9AC9 | na | | Relay Checksum | base 16 # | na | 9AC9 | na |
| | Model Type | coded # | tyP | 4 | 4, 20, 30, 5, 55, 66 | | Model Type | coded # | tyP | 5 | 4, 20, 30, 5, 55, 66 |
| | Configuration | coded # | CFg | 01 | 00, 01 | | Configuration | coded # | CFg | 01 | 00, 01 |
| | Cubic Foot | cu ft | CuFt | (cab size) | 12, 23, 30, 50 | | Cubic Foot | cu ft | CuFt | (cab size) | 12, 23, 30, 45, 50 |
| | Pull Down Offset | °C | Pd oFSt | 0.0 | +/- 5.0 | | Pull Down Offset | °C | Pd oFSt | 0.5 | +/- 5.0 |
| | Control Offset | °C | Cnt oFSt | 0.2 | +/- 5.0 | | Control Offset | °C | Cnt oFSt | 0.0 | +/- 5.0 |
| | Display Offset | °C | oFSt | -1.4 | +/- 3.0 | | Display Offset | °C | oFSt | -0.9 | +/- 3.0 |
| | Control UCL | °C | Cnt ucl | 0.8 | 0 to 5 | | Control UCL | °C | Cnt ucl | 0.8 | 0 to 5 |
| | Control LCL | °C | Cnt lcl | -0.2 | -5 to 0 | | Control LCL | °C | Cnt lcl | -0.2 | -5 to 0 |
| | Defrost Interval | hours | dEF Int | 1 | 1 to 12 | | Defrost Interval | hours | dEF Int | 1 | 1 to 12 |
| | Defrost Duration | minutes | dEF dur | 15 | 15 to 30 | | Defrost Duration | minutes | dEF dur | 15 | 15 to 30 |
| | Defr. Hi (temp out) | °C | dEF Hi | 5 | 1 to 15 | | Defr. Hi (temp out) | °C | dEF Hi | 5 | 1 to 15 |
| | Quality Mode | on/off | qUA | 00 | 00, 01 | | Quality Mode | on/off | qUA | 00 | 00, 01 |
| | BOT Request | on/off | bot | 00 | 00, 01 | | BOT Request | on/off | bot | 00 | 00, 01 |
| | Perimeter Heater PWM | % | Per | 35 | 20 to 100 | | Perimeter Heater PWM | % | Per | 35 | 20 to 100 |
| | Cntr.Probe | °C | Cnt Prb | °C | reading | | Cntr.Probe | °C | Cnt Prb | °C | reading |
| | Amb. Probe | °C | Anb Prb | °C | reading | | Amb. Probe | °C | Anb Prb | °C | reading |
| | Lower Display Probe | °C | d lo Prb | °C | reading | | Lower Display Probe | °C | d lo Prb | °C | reading |
| | Defr.Prb | °C | dEF Prb | °C | reading | | Defr.Prb | °C | dEF Prb | °C | reading |

Suffix 17 G54 Rev. C Paramters Apr. 2019 - Feb. 2020

| "-20 Freezer" | | | | | | "-30 Freezer" | | | | | |
|---------------|----------------------|-----------|----------------|------------|---------------|---------------|----------------------|-----------|----------------|------------|----------------------|
| 20 | System Parameter | Units | Display "text" | Value | Range | 30 | System Parameter | Units | Display "text" | Value | Range |
| PROG | Set point | °C | na | -20.0 | -25 to -15 | PROG | Set point | °C | na | -30.0 | -35 to -15 |
| MODE | Warm Alarm | °C | na | -10.0 | (sp+5) to -10 | MODE | Warm Alarm | °C | na | -20.0 | (sp+5) to -10 |
| | Cold Alarm | °C | na | -25.0 | -35 to (sp-5) | | Cold Alarm | °C | na | -35.0 | -40 to (sp-5) |
| | "SER" | na | SEr | na | na | | "SER" | na | SEr | na | na |
| SERVICE | Temp of Unit | °C | na | 20 | 4,20,30,5 | SERVICE | Temp of Unit | °C | na | 30 | 4,20,30,5 |
| MODE | CPU Checksum | base 16 # | na | A678 | na | MODE | CPU Checksum | base 16 # | na | A678 | na |
| | Relay Checksum | base 16 # | na | 9AC9 | na | | Relay Checksum | base 16 # | na | 9AC9 | na |
| | Model Type | coded # | tyP | 20 | 20, 30 | | Model Type | coded # | tyP | 30 | 4, 20, 30, 5, 55, 66 |
| | Configuration | coded # | CFg | 02 | 00, 01, 02 | | Configuration | coded # | CFg | 01 | 00, 01 |
| | Cubic Foot | cu ft | CuFt | (cab size) | 23, 30 | | Cubic Foot | cu ft | CuFt | (cab size) | 12, 23, 30, 50 |
| | Pull Down Offset | °C | Pd oFSt | -7.0 | +/- 10 | | Pull Down Offset | °C | Pd oFSt | -2.0 | +/- 5.0 |
| | Control Offset | °C | Cnt oFSt | -3.0 | +/- 10 | | Control Offset | °C | Cnt oFSt | -0.5 | +/- 5.0 |
| | Display Offset | °C | oFSt | -2.2 | +/- 10 | | Display Offset | °C | oFSt | 1.3 | +/- 10 |
| | Control UCL | °C | Cnt ucl | 0.5 | 0 to 5 | | Control UCL | °C | Cnt ucl | 0.5 | 0 to 5 |
| | Control LCL | °C | Cnt lcl | -0.5 | -5 to 0 | | Control LCL | °C | Cnt lcl | -0.5 | -5 to 0 |
| | Quality Mode | on/off | qUA | 00 | 00, 01 | | Adaptive Defrost | on/off | AdP dEF | 01 | 00, 01 |
| | BOT Request | on/off | bot | 00 | 00, 01 | | Defrost Interval | hours | dEF Int | 6 | 4 to 8 |
| | Perimeter Heater PWM | % | Per | 65 | 20 to 100 | | Defrost Duration | minutes | dEF dur | 20 | 20 to 30 |
| | Cntr.Probe | °C | Cnt Prb | °C | reading | | Defr. Hi (temp out) | °C | dEF Hi | 1 | 1 to 10 |
| | Amb. Probe | °C | Anb Prb | °C | reading | | Drip Pan temp out | °C | Pan Hi | 3 | 0 to 10 |
| | | | | | | | Quality Mode | on/off | qUA | 00 | 00, 01 |
| | | | | | | | BOT Request | on/off | bot | 00 | 00, 01 |
| | | | | | | | Perimeter Heater PWM | % | Per | 90 | 20 to 100 |
| | | | | | | | Cntr.Probe | °C | Cnt Prb | °C | reading |
| | | | | | | | Amb. Probe | °C | Anb Prb | °C | reading |
| | | | | | | | Lower Display Probe | °C | d lo Prb | °C | reading |
| | | | | | | | Defr.Prb | °C | dEF Prb | °C | reading |
| | | | | | | | Drip Pan Probe | °C | PAn Prb | °C | reading |

Suffix 17, 18, 19 G54 Rev. D Paramters Feb. 2020 - May 2021

| "+4 Refrigerator" | | | | | | "+5 Refrigerator" | | | | | |
|-------------------|----------------------|-----------|----------------|-------------|----------------------|-------------------|----------------------|-----------|----------------|-------------|----------------------|
| 4 | System Parameter | Units | Display "text" | Value | Range | 5 | System Parameter | Units | Display "text" | Value | Range |
| PROG | Set point | °C | na | 4.0 | 3 to 7 | PROG | Set point | °C | na | 5.0 | 3 to 7 |
| MODE | Warm Alarm | °C | na | 5.5 | (sp+1) to 15 | MODE | Warm Alarm | °C | na | 8.0 | (sp+1) to 15 |
| | Cold Alarm | °C | na | 2.0 | 0 to (sp-1) | | Cold Alarm | °C | na | 2.0 | 0 to (sp-1) |
| | "SER" | na | SEr | na | na | | "SER" | na | SEr | na | na |
| SERVICE | Temp of Unit | °C | na | 4 | 4,20,30,5 | SERVICE | Temp of Unit | °C | na | 5 | 4,20,30,5 |
| MODE | CPU Checksum | base 16 # | na | A901 | na | MODE | CPU Checksum | base 16 # | na | A901 | na |
| | Relay Checksum | base 16 # | na | D284 / 3429 | na | | Relay Checksum | base 16 # | na | D284 / 3429 | na |
| | Model Type | coded # | tyP | 4 | 4, 20, 30, 5, 55, 66 | | Model Type | coded # | tyP | 5 | 4, 20, 30, 5, 55, 66 |
| | Configuration | coded # | CFg | 01 | 00, 01 | | Configuration | coded # | CFg | 01 | 00, 01 |
| | Cubic Foot | cu ft | CuFt | (cab size) | 12, 23, 30, 50 | | Cubic Foot | cu ft | CuFt | (cab size) | 12, 23, 30, 45, 50 |
| | Pull Down Offset | °C | Pd oFSt | 0.0 | +/- 5.0 | | Pull Down Offset | °C | Pd oFSt | 0.5 | +/- 5.0 |
| | Control Offset | °C | Cnt oFSt | 0.2 | +/- 5.0 | | Control Offset | °C | Cnt oFSt | 0.0 | +/- 5.0 |
| | Display Offset | °C | oFSt | -1.4 | +/- 3.0 | | Display Offset | °C | oFSt | -0.9 | +/- 3.0 |
| | Control UCL | °C | Cnt ucl | 0.8 | 0 to 5 | | Control UCL | °C | Cnt ucl | 0.8 | 0 to 5 |
| | Control LCL | °C | Cnt lcl | -0.2 | -5 to 0 | | Control LCL | °C | Cnt lcl | -0.2 | -5 to 0 |
| | Door Ajar | minutes | door ajar | 3 | 1, 2, 3 | | Door Ajar | minutes | door ajar | 3 | 1, 2, 3 |
| | Defrost Interval | hours | dEF Int | 1 | 1 to 12 | | Defrost Interval | hours | dEF Int | 1 | 1 to 12 |
| | Defrost Duration | minutes | dEF dur | 15 | 15 to 30 | | Defrost Duration | minutes | dEF dur | 15 | 15 to 30 |
| | Defr. Hi (temp out) | °C | dEF Hi | 5 | 1 to 15 | | Defr. Hi (temp out) | °C | dEF Hi | 5 | 1 to 15 |
| | Quality Mode | on/off | qUA | 00 | 00, 01 | | Quality Mode | on/off | qUA | 00 | 00, 01 |
| | BOT Request | on/off | bot | 00 | 00, 01 | | BOT Request | on/off | bot | 00 | 00, 01 |
| | Perimeter Heater PWM | % | Per | 35 | 20 to 100 | | Perimeter Heater PWM | % | Per | 35 | 20 to 100 |
| | Cntr.Probe | °C | Cnt Prb | °C | reading | | Cntr.Probe | °C | Cnt Prb | °C | reading |
| | Amb. Probe | °C | Anb Prb | °C | reading | | Amb. Probe | °C | Anb Prb | °C | reading |
| | Lower Display Probe | °C | d lo Prb | °C | reading | | Lower Display Probe | °C | d lo Prb | °C | reading |
| | Defr.Prb | °C | dEF Prb | °C | reading | | Defr.Prb | °C | dEF Prb | °C | reading |

| Suffix 17, 18, 19 G54 Rev. D Paramters Feb. 2020 - May 2021 | | | | | | | | | | | |
|---|----------------------|-----------|----------------|-------------|----------------|---------------|----------------------|-----------|----------------|-------------|----------------------|
| "-20 Freezer" | | | | | | "-30 Freezer" | | | | | |
| 20 | System Parameter | Units | Display "text" | Value | Range | 30 | System Parameter | Units | Display "text" | Value | Range |
| PROG | Set point | °C | na | -20.0 | -25 to -15 | PROG | Set point | °C | na | -30.0 | -35 to -15 |
| MODE | Warm Alarm | °C | na | -10.0 | (sp+5) to -10 | MODE | Warm Alarm | °C | na | -20.0 | (sp+5) to -10 |
| | Cold Alarm | °C | na | -25.0 | -35 to (sp-5) | | Cold Alarm | °C | na | -35.0 | -40 to (sp-5) |
| | "SER" | na | SEr | na | na | | "SER" | na | SEr | na | na |
| SERVICE | Temp of Unit | °C | na | 20 | 4,20,30,5 | SERVICE | Temp of Unit | °C | na | 30 | 4,20,30,5 |
| MODE | CPU Checksum | base 16 # | na | A901 | na | MODE | CPU Checksum | base 16 # | na | A901 | na |
| | Relay Checksum | base 16 # | na | D284 / 3429 | na | | Relay Checksum | base 16 # | na | D284 / 3429 | na |
| | Model Type | coded # | tyP | 20 | 20, 30 | | Model Type | coded # | tyP | 30 | 4, 20, 30, 5, 55, 66 |
| | Configuration | coded # | CFg | 03 | 00, 01, 02, 03 | | Configuration | coded # | CFg | 01 | 00, 01 |
| | Cubic Foot | cu ft | CuFt | (cab size) | 23, 30 | | Cubic Foot | cu ft | CuFt | (cab size) | 12, 23, 30, 50 |
| | Pull Down Offset | °C | Pd oFSt | -2.0 | +/- 10 | | Pull Down Offset | °C | Pd oFSt | -2.0 | +/- 5.0 |
| | Control Offset | °C | Cnt oFSt | -0.8 | +/- 10 | | Control Offset | °C | Cnt oFSt | -0.5 | +/- 5.0 |
| | Display Offset | °C | oFSt | -2.2 | +/- 10 | | Display Offset | °C | oFSt | 1.3 | +/- 10 |
| | Control UCL | °C | Cnt ucl | 0.1 | 0 to 5 | | Control UCL | °C | Cnt ucl | 0.5 | 0 to 5 |
| | Control LCL | °C | Cnt lcl | -0.1 | -5 to 0 | | Control LCL | °C | Cnt lcl | -0.5 | -5 to 0 |
| | Door Ajar | minutes | door ajar | 03 | 01, 02, 03 | | Door Ajar | minutes | door ajar | 3 | 1, 2, 3 |
| | Quality Mode | on/off | qUA | 00 | 00, 01 | | Adaptive Defrost | on/off | AdP dEF | 01 | 00, 01 |
| | BOT Request | on/off | bot | 00 | 00, 01 | | Defrost Interval | hours | dEF Int | 6 | 4 to 8 |
| | Perimeter Heater PWM | % | Per | 65 | 20 to 100 | | Defrost Duration | minutes | dEF dur | 20 | 20 to 30 or na |
| | Cntr.Probe | °C | Cnt Prb | °C | reading | | Defr. Hi (temp out) | °C | dEF Hi | 1 | 1 to 10 |
| | Amb. Probe | °C | Anb Prb | °C | reading | | Drip Pan temp out | °C | Pan Hi | 3 | 0 to 10 |
| | | | | | | | Quality Mode | on/off | qUA | 00 | 00, 01 |
| | | | | | | | BOT Request | on/off | bot | 00 | 00, 01 |
| | | | | | | | Perimeter Heater PWM | % | Per | 90 | 20 to 100 |
| | | | | | | | Cntr.Probe | °C | Cnt Prb | °C | reading |
| | | | | | | | Amb. Probe | °C | Anb Prb | °C | reading |
| | | | | | | | Lower Display Probe | °C | d lo Prb | °C | reading |
| | | | | | | | Defr.Prb | °C | dEF Prb | °C | reading |
| | | | | | | | Drip Pan Probe | °C | PAn Prb | °C | reading |

Suffix 17, 18, 19 G54 Rev. E Paramters May 2021 to Current

| "+4 Refrigerator" | | | | | | "+5 Refrigerator" | | | | | |
|-------------------|----------------------|-----------|----------------|------------|----------------------|-------------------|----------------------|-----------|----------------|------------|----------------------|
| 4 | System Parameter | Units | Display "text" | Value | Range | 5 | System Parameter | Units | Display "text" | Value | Range |
| PROG | Set point | °C | na | 4.0 | 3 to 7 | PROG | Set point | °C | na | 5.0 | 3 to 7 |
| MODE | Warm Alarm | °C | na | 5.5 | (sp+1) to 15 | MODE | Warm Alarm | °C | na | 8.0 | (sp+1) to 15 |
| | Cold Alarm | °C | na | 2.0 | 0 to (sp-1) | | Cold Alarm | °C | na | 2.0 | 0 to (sp-1) |
| | "SER" | na | SEr | na | na | | "SER" | na | SEr | na | na |
| SERVICE | Temp of Unit | °C | na | 4 | 4,20,30,5 | SERVICE | Temp of Unit | °C | na | 5 | 4,20,30,5 |
| MODE | CPU Checksum | base 16 # | na | 36DE | na | MODE | CPU Checksum | base 16 # | na | 36DE | na |
| | Relay Checksum | base 16 # | na | 9862 | na | | Relay Checksum | base 16 # | na | 9862 | na |
| | Model Type | coded # | tyP | 4 | 4, 20, 30, 5, 55, 66 | | Model Type | coded # | tyP | 5 | 4, 20, 30, 5, 55, 66 |
| | Configuration | coded # | CFg | 04 | 00, 01 | | Configuration | coded # | CFg | 04 | 00, 01 |
| | Cubic Foot | cu ft | CuFt | (cab size) | 12, 23, 30, 50 | | Cubic Foot | cu ft | CuFt | (cab size) | 12, 23, 30, 45, 50 |
| | Pull Down Offset | °C | Pd oFSt | 0.0 | +/- 5.0 | | Pull Down Offset | °C | Pd oFSt | 0.5 | +/- 5.0 |
| | Control Offset | °C | Cnt oFSt | 0.2 | +/- 5.0 | | Control Offset | °C | Cnt oFSt | 0.0 | +/- 5.0 |
| | Display Offset | °C | oFSt | -1.4 | +/- 3.0 | | Display Offset | °C | oFSt | -0.9 | +/- 3.0 |
| | Control UCL | °C | Cnt ucl | 0.8 | 0 to 5 | | Control UCL | °C | Cnt ucl | 0.8 | 0 to 5 |
| | Control LCL | °C | Cnt lcl | -0.2 | -5 to 0 | | Control LCL | °C | Cnt lcl | -0.2 | -5 to 0 |
| | Door Ajar | minutes | door ajar | 3 | 1, 2, 3 | | Door Ajar | minutes | door ajar | 3 | 1, 2, 3 |
| | Defrost Interval | hours | dEF Int | 1 | 1 to 12 | | Defrost Interval | hours | dEF Int | 1 | 1 to 12 |
| | Defrost Duration | minutes | dEF dur | 15 | 15 to 30 | | Defrost Duration | minutes | dEF dur | 15 | 15 to 30 |
| | Defr. Hi (temp out) | °C | dEF Hi | 5 | 1 to 15 | | Defr. Hi (temp out) | °C | dEF Hi | 5 | 1 to 15 |
| | Quality Mode | on/off | qUA | 00 | 00, 01 | | Quality Mode | on/off | qUA | 00 | 00, 01 |
| | BOT Request | on/off | bot | 00 | 00, 01 | | BOT Request | on/off | bot | 00 | 00, 01 |
| | Perimeter Heater PWM | % | Per | 35 | 20 to 100 | | Perimeter Heater PWM | % | Per | 35 | 20 to 100 |
| | Cntr.Probe | °C | Cnt Prb | °C | reading | | Cntr.Probe | °C | Cnt Prb | °C | reading |
| | Amb. Probe | °C | Anb Prb | °C | reading | | Amb. Probe | °C | Anb Prb | °C | reading |
| | Lower Display Probe | °C | d lo Prb | °C | reading | | Lower Display Probe | °C | d lo Prb | °C | reading |
| | Defr.Prb | °C | dEF Prb | °C | reading | | Defr.Prb | °C | dEF Prb | °C | reading |

Suffix 17, 18, 19 G54 Rev. E Paramters May 2021 to Current

| "-20 Freezer" | | | | | | "-30 Freezer" | | | | | |
|---------------|----------------------|-----------|----------------|------------|--------------------|---------------|----------------------|-----------|----------------|------------|----------------------|
| 20 | System Parameter | Units | Display "text" | Value | Range | 30 | System Parameter | Units | Display "text" | Value | Range |
| PROG | Set point | °C | na | -20.0 | -25 to -15 | PROG | Set point | °C | na | -30.0 | -35 to -15 |
| MODE | Warm Alarm | °C | na | -10.0 | (sp+5) to -10 | MODE | Warm Alarm | °C | na | -20.0 | (sp+5) to -10 |
| | Cold Alarm | °C | na | -25.0 | -35 to (sp-5) | | Cold Alarm | °C | na | -35.0 | -40 to (sp-5) |
| | "SER" | na | SEr | na | na | | "SER" | na | SEr | na | na |
| SERVICE | Temp of Unit | °C | na | 20 | 4,20,30,5 | SERVICE | Temp of Unit | °C | na | 30 | 4,20,30,5 |
| MODE | CPU Checksum | base 16 # | na | 36DE | na | MODE | CPU Checksum | base 16 # | na | 36DE | na |
| | Relay Checksum | base 16 # | na | 9862 | na | | Relay Checksum | base 16 # | na | 9862 | na |
| Starts here | Model Type | coded # | tyP | 20 | 20, 30 | | Model Type | coded # | tyP | 30 | 4, 20, 30, 5, 55, 66 |
| | Configuration | coded # | CFg | 04 | 00, 01, 02, 03, 04 | | Configuration | coded # | CFg | 04 | 00, 01 |
| | Cubic Foot | cu ft | CuFt | (cab size) | 23, 30 | | Cubic Foot | cu ft | CuFt | (cab size) | 12, 23, 30, 50 |
| | Pull Down Offset | °C | Pd oFSt | -2.0 | +/- 10 | | Pull Down Offset | °C | Pd oFSt | -2.0 | +/- 5.0 |
| | Control Offset | °C | Cnt oFSt | -0.8 | +/- 10 | | Control Offset | °C | Cnt oFSt | -0.5 | +/- 5.0 |
| | Display Offset | °C | oFSt | -2.2 | +/- 10 | | Display Offset | °C | oFSt | 1.3 | +/- 10 |
| | Control UCL | °C | Cnt ucl | 0.1 | 0 to 5 | | Control UCL | °C | Cnt ucl | 0.5 | 0 to 5 |
| | Control LCL | °C | Cnt lcl | -0.1 | -5 to 0 | | Control LCL | °C | Cnt lcl | -0.5 | -5 to 0 |
| | Door Ajar | minutes | door ajar | 03 | 01, 02, 03 | | Door Ajar | minutes | door ajar | 3 | 1, 2, 3 |
| | Quality Mode | on/off | qUA | 00 | 00, 01 | | Adaptive Defrost | on/off | AdP dEF | 01 | 00, 01 |
| | BOT Request | on/off | bot | 00 | 00, 01 | | Defrost Interval | hours | dEF Int | 6 | 4 to 8 |
| | Perimeter Heater PWM | % | Per | 65 | 20 to 100 | | Defrost Duration | minutes | dEF dur | 20 | 20 to 30 or na |
| | Cntr.Probe | °C | Cnt Prb | °C | reading | | Defr. Hi (temp out) | °C | dEF Hi | 1 | 1 to 10 |
| | Amb. Probe | °C | Anb Prb | °C | reading | | Drip Pan temp out | °C | Pan Hi | 3 | 0 to 10 |
| | | | | | | | Quality Mode | on/off | qUA | 00 | 00, 01 |
| | | | | | | | BOT Request | on/off | bot | 00 | 00, 01 |
| | | | | | | | Perimeter Heater PWM | % | Per | 90 | 20 to 100 |
| | | | | | | | Cntr.Probe | °C | Cnt Prb | °C | reading |
| | | | | | | | Amb. Probe | °C | Anb Prb | °C | reading |
| | | | | | | | Lower Display Probe | °C | d lo Prb | °C | reading |
| | | | | | | | Defr.Prb | °C | dEF Prb | °C | reading |
| | | | | | | | Drip Pan Probe | °C | PAn Prb | °C | reading |

TSX Sensor Matrix (RoHS compliant)

| Part Number | Description | Location | Sensor Color |
|-------------|----------------------------------|--|--------------|
| 309159H09 | Display Sensor, (PT1000 Ohm RTD) | Refrigerator and Freezer Bottle | Green |
| 309159H10 | Display Sensor, (PT1000 Ohm RTD) | -20°C Freezer Bottle | Green |
| 312096H06 | Display Sensor, (LM335) | Blood Bank +4°C Lower Bottle | Yellow |
| 326737H05 | Defrost Sensor, (PT1000 Ohm RTD) | Refrigerator and Freezer Evaporator | Blue |
| 326737H03 | Control Sensor, (PT1000 Ohm RTD) | Refrigerator and Freezer Evaporator On Wall of Early -20°C Freezers | Red |
| 326737H04 | Control Sensor, (PT1000 Ohm RTD) | On Wall of -20°C Freezers | Red |
| 326739H01 | Ambient Sensor, (LM335) | Refrigerator and Freezer Electrical Deck | White |
| 326739H03 | Drip Pan Sensor, (LM335) | Freezer Evaporator Pan | Black |
| 85089G10 | 100 Ohm RTD | 4-20mA and Chart Recorder Sensor | Grey |


Display Sensor and Controller Accuracy

| Model | Component | Part Number | Description | Temperature Range | Accuracy Specification | Resolution |
|------------|---------------------------|-------------|----------------------|--------------------------|---|----------------|
| TSX1230 | Display Temperature Probe | 309159h09 | Platinum RTD / 1000Ω | -150°C to 60°C | +/-0.30° C at 0° C, Tested to DIN EN 60751, Class B | Not Applicable |
| TSX2330 | | | | | | |
| TSX3030 | | | | | | |
| TSX5030 | | 309159h10 | | | | |
| TSX2320 | | | | | | |
| TSX3020 | | | | | | |
| TSX1204 | | 309159h09 | | | | |
| TSX2304 | | | | | | |
| TSX3004 | | | | | | |
| TSX5004 | | 309159h10 | | | | |
| TSX1205 | | 309159h09 | | | | |
| TSX2305 | | | | | | |
| TSX3005 | | | | | | |
| TSX4505 | | 309159h10 | | | | |
| TSX5005 | | | | | | |
| All Models | Controller | 304864H06 | Microprocessor | +99.9°C to -99.9°C | ± 0.5°C | 0.1°C |


| Model | Component | Part Number | Description | Temperature Range | Accuracy Specification | Resolution |
|------------|---------------------------|-------------|----------------------|--------------------------|---|----------------|
| TSX1230 | Temperature Control Probe | 326737H03 | Platinum RTD / 1000Ω | -50°C to 300°C | +/-0.15° C at 0° C, Tested to DIN EN 60751, Class A | Not Applicable |
| TSX2330 | | | | | | |
| TSX3030 | | | | | | |
| TSX5030 | | | | | | |
| TSX2320 | | 326737H04 | | | | |
| TSX3020 | | | | | | |
| TSX1204 | | 326737H03 | | | | |
| TSX2304 | | | | | | |
| TSX3004 | | | | | | |
| TSX5004 | | | | | | |
| TSX1205 | | | | | | |
| TSX2305 | | | | | | |
| TSX3005 | | | | | | |
| TSX4505 | | | | | | |
| TSX5005 | | | | | | |
| All Models | Controller | 304864H06 | Microprocessor | +99.9°C to -99.9°C | ± 0.5°C | 0.1°C |

Standard CPU Sensor / Harness Connections


+4°C Blood Bank Refrigerator




+5°C Pharmacy and General-Purpose Refrigerator





-30°C Freezer – CFg 00 Electric Defrost




-30°C Freezer – CFg 01, 04 Hot Gas Defrost










Sensor Locations Inside of Various TSX Units

Display Bottle: (green) all refrigerators and freezers except -20 cold wall, top left side wall.


Note: *Sensor bottle is now shipped with a solid cap to reduce spillage. The sensor must be installed by the customer or service tech.*





BB Display and Lower Bottle: (Green and Yellow) Blood banks only – lower right on the back wall. Also Dixell overcold protection sensor in the bottle with the Yellow sensor.



Refrigerator Control: (red) “L” bracket inside the evaporator in front of the coil.



Defrost: (blue) clamped to the evaporator outlet tube (varies with the different refrigerator / freezer evaporator sizes), or in high door opening cases, in the fins of the coil.



-20 Control: (red) mounted to the left side wall with a metal bracket. The bracket must be thermally isolated from the wall by nylon spacers.

***Note:** in this picture, the normal spacers have been replaced by longer spacers to clarify their location.*




-20 Control Config = 03-04 Suffix 19 and later: (red) mounted to the Lower left side wall with a plastic bracket. The bracket must be thermally isolated from the wall by nylon spacers. Sensor mounting location is lower on 3020's than the 2320's due to the extra Evaporator loop.



Control sensor mounted inside the sensor cover.

CAD view of evap lines and sensor



PEP


Display Sensor

Note: This will not work, sensor(s) must be inserted into the bottle.



Inverter Information

Embraco High Performance (HP)




| Inverter Technical Specifications | |
|--|--|
| Nominal input voltage range | 220 VAC to 240 VAC (220 / 240 V model) |
| Operating input voltage range | 150 VAC to 264 VAC (220 / 240 V model) |
| AC input frequency | 50 Hz to 60 Hz |
| Speed range | 2000 to 4500 rpm |
| Limits of ambient temperature | +43°C to -20°C |
| Maximum storage external relative humidity | 85% |
| Air forced ventilation (min) | 2 m/s |


Embraco High Performance (HP) LED Codes

| LED Status | LED Status Description |
|-------------------------------|--|
| 1 flash 15 seconds period | Normal operation |
| 1 flash 0.5 seconds period | Over temperature limit |
| 2 flashes 5 seconds period | Check control signal |
| 3 flashes 5 seconds period | Check inverter |
| 4 flashes 5 seconds period | Check compressor |
| 5 flashes 5 seconds period | Input voltage (under 150 VAC or over 264 VAC) |

Maia Inverter (Started with some Suffix 17 Models)



INVERTER COVER NOT
SHOWN IN THIS VIEW



No Transformer on 2304 / 5 Refrigerators.

Maia LED Codes

| Maia Inverter Troubleshooting Flash Codes | | | |
|---|------------|--------|--|
| LED Status | Period | Color | Description |
| 1 Flash | 30 seconds | Green | Normal Operation |
| 2 Flashes | 5 seconds | Green | Communication Problem |
| 3 Flashes | 5 seconds | Red | Inverter Problem |
| 4 Flashes | 5 seconds | Orange | Compressor Problem |
| No Flash | - | - | No input power / Damaged inverter / Dixell cut off input power due to anti-freeze protection |


MAIA flash sequences – normal is now 30 sec as opposed to the HP 15 sec

Condenser Fan

- This unit can operate the condenser fan at either Low Speed, or High-Speed.
- The fan speed is determined based on three conditions below:
 - Ambient temperature (Greater than 32°C sets the fan to high speed)
 - Temperature at ambient probe (Greater than 32°C for freezers, and 40°C for refrigerators sets the fan to High-Speed.
- The compressor speed (Greater than, or equal to 4600RPM will result in High-Speed fan operation.)

| Current Speed | Scenario | New Speed |
|---------------|---------------------------------------|------------|
| Low Speed | Ambient Temp \geq Threshold | High Speed |
| High Speed | Ambient Temp < Threshold - Hysteresis | Low Speed |

Note: If the Ambient Sensor fails the condenser fan will go to the High-Speed based on thresholds of 40°C for refrigerators, and 32°C for -20°, -30°C freezers.



K1 = Line voltage, K3 = Low speed, n1=2100 rpm, n2 = 1500 rpm





Electric Condensate Heater 328840H0X

Used on 12', 23', 30', 45' and 50' TSX units since 2019.



To test - check the amp draw without water, and then with water. It will start dry around .3 Amps, and will climb with water to around 1.1 amps.



Chromatography with Safety Outlet

Modified from the User Manual –

Chromatography Refrigerators are equipped with a convenience outlet located on the back wall. Power to the outlet is protected by a 5A circuit breaker. If there is a fault, the circuit can be reset by toggling the Chromatography Power Inlet Switch (Item B in *Figure 2*). **The total load of the equipment using the convenience outlet shall not exceed 4 Amps.**

WARNING: The convenience outlet may be *LIVE* when the unit is connected to an electrical power source, regardless of the position of the key switch.

The chromatography refrigerator outlet is also equipped with a safety circuit. In the event the safety circuit is tripped (the safety sensor senses a flammable gas), the power to the internal outlet will be discontinued and the user will need to manually reset the safety circuit assembly using the reset switch at the back of the unit before power is restored. To do this, toggle the switch labeled “Switch Chromatography Manual Reset” off and then on. There will be a power delay to the outlet for approximately 4 minutes where the red LED will be illuminated. When power is available at the outlet, the red LED will turn off and the green LED will illuminate. If the circuit trips from RED back to GREEN repeatedly or never turns GREEN, replace the two circuit boards

Flammable materials should not be stored in the refrigerator.

The safety circuit has two LEDs:

- **RED** (on the top) - When the red LED is ON, it indicates there is no power to the AC outlet inside the unit.
- **GREEN** (on the bottom) - When the green LED is ON, it indicates there is power to the AC outlet inside the unit.



For the US, the NEMA 5-15 GFCI Receptacle has a standard self-test feature with color indicators for status:

- **Solid Red** - Indicates the breaker has opened the circuit. If this occurs, unplug the equipment, check the circuit, and reset the circuit breaker.
- **Flashing Red** - Indicates “End of Life” and that the receptacle needs to be replaced.
- **Solid Green** - Indicates the outlet is functional.



NEMA 5-15 Receptacle Lights


The first-time power is applied to a chromatography unit with the NEMA 15 outlet installed as opposed to the European, British, or Chinese outlets, it will need the RESET button pressed in before the outlet will power up. This step should not be needed during sensor trips and switch position changes at the back of the unit.

The chromatography refrigerator has two power inlets at the back of the unit:

1. **Main Power Inlet (A)** - This is used to connect the AC mains power cord to the chromatography refrigerator.
2. **Chromatography Power Inlet (C)** - This is used to connect the AC mains power cord to the chromatography outlet assembly inside of the cabinet.

The chromatography refrigerator has three switches at the back of the unit:

1. **Mains Power Inlet Switch (B)** - This is used to turn ON or OFF the AC supply to the refrigerator.
2. **Chromatography Power Inlet Switch (D)** - This is used to turn ON or OFF the AC inlet supply to the receptacle inside the refrigerator.
3. **Chromatography Manual Reset Switch (E)** - This is used to manually reset the safety circuit in the AC receptacle assembly after the LED turns Red from Green due to a safety trip event which disengages power to the outlet inside the refrigerator.



Chromatography Refrigerator Power Inlet and Switches

Refer to the label to the right of the switches at the back of the refrigerator for identification.

Due to the nature of new HC leak detecting sensor, there are additional considerations when cleaning a chromatography unit.

The refrigerator is designed to be cleaned using a mild detergent (such as Lysol® wipes (non-aerosol) or Formula 409®) and water. Lightly spray the interior storage components and wipe them dry with a soft cloth or spray the cloth first and then wipe interior surfaces. Do not spray directly on the center outlet column. Use a damp cloth to clean around the outlet. Use of other cleaners may contain chemicals that will turn off power to the outlet. If this occurs, continue to clean the refrigerator, and leave the doors open for approximately 5 minutes to allow any vapors to exit the refrigerator.

If the outlet power turns off due to chemicals, the red light will illuminate, and the user will need to manually reset the safety circuit using the “Switch Manual Reset Chromatography” at the back of the unit before power is restored. Sufficient time will be needed to dissipate the chemical vapors that tripped the sensor. The door can be opened to help dissipate the chemical vapors or an external fan can be used if the sensor is being tested before the unit is placed into normal operation. Toggling this switch off then on will restart the safety system and the green light should illuminate after approximately 4 minutes. If the light repeatedly changes from red to green or the red-light stays illuminated, ensure there are no chemical vapors present in the refrigerator and reset the system. If this continues, replace the boards in the outlet circuit.

The sensor can be tested with small amounts of R290 propane, but is also sensitive to iso-butane and LPG. Other gases that may trip the sensor are LNG, alcohol and smoke.

Other reasons for replacing the outlet sensor:

- 1) The sensor needs replacement if the temperature goes below freezing in the cabinet. Use system alarms to ensure the temperature inside the unit is always above 0°C. If exposed to freezing conditions the sensor shall be replaced.
- 2) Simply time – **the sensor shall be replaced every five years**. This is due to life expectancy regardless of any exposure to freezing, unintentional or valid trips.


If the convenience outlet has no power :


Check for all switches in the “ON” position (“1”) on the rear of the unit (pg 75). Try cycling the Chromatography Manual Reset Switch and the Chromatography Power Inlet Switch to OFF and then ON. For units with the GFCI receptacle, make sure the RESET switch is pressed in.

If the sensor is tripped and the light is red, check for chemicals that the customer may have stored in the refrigerator in an unsealed container (ethanol or alcohol based). Leave the doors open for approximately 5 minutes to allow any vapors to exit the refrigerator and cycle the Chromatography Manual Reset Switch off and back on. Again, 4 minute delay before the light will go green.

If the lights repeatedly change from red to green or the red lights stays illuminated after trying to reset, replace the sensor board before replacing the control board.

HC Sensor PPM Sensativity Graph





Blood Bank Anti-Freeze Protection

1 Objective

Dixell XR02CX is installed with compressor relay contact in series with Relay Board compressor output. The Dixell is configured such that when the control probe temperature drops below 2.8C the compressor is interrupted until the control probe temperature recovers to above 4.0C. The evaporator probe and defrost function is disabled.

2 Dixell XR02CX Settings

The temperature set point is displayed by briefly pressing the SET button – record this value The

Probe Fail action is set up to stop the compressor on probe failure.

To access other parameters press the SET button and button until HY is displayed.

Release the buttons then press the SET button and button again until Pr2 is displayed.

The first parameter is HY; press SET to display its value then press SET again to move to the next parameter. Repeat until all parameters have been displayed and recorded.

To access the set point press the SET button until the units symbol “C” flashes then use the up and down buttons to change the value.

Set the parameters in the order listed in the table, changing the set point last.

| Order | Parameter description | Label | Set | Meaning |
|-------|--|------------|------|-------------------------------|
| 27 | Set point | Set | 2.8 | Output turns off below 2.8C |
| 17 | Differential | HY | 1.2 | Output turns on above 4.0C |
| 18 | Minimum set point | LS | -4.0 | |
| 19 | Maximum set point | US | 9.9 | |
| 20 | Control probe calibration | ot | 0.0 | |
| 21 | Evaporator probe presence | P2 | n | Probe not present |
| 22 | Evaporator probe calibration | oE | 0.0 | |
| 23 | Outputs delay at start up | od | 0 | |
| 24 | Anti-short cycle delay | AC | 0 | |
| 25 | Compressor ON with faulty probe | CY | 0 | USE THIS FOR PROBE FAIL = OFF |
| 26 | Compressor OFF time with faulty probe | Cn | 30 | USE THIS FOR PROBE FAIL = OFF |
| 1 | Temperature measurement unit | CF | C | |
| 2 | Resolution | rE | dE | |
| 3 | Probe displayed | Ld | P1 | Control Probe |
| 4 | Display temperature delay | dY | 0 | Display is immediate no delay |
| 5 | Defrost termination temperature | dE | 4.6 | |
| 6 | Interval between defrost cycles | id | 0 | 0 = Defrost disabled |
| 7 | Maximum length for defrost | Nd | 0 | 0 = Defrost disabled |
| 8 | Displaying during defrost | dF | rt | Control probe |
| 9 | Maximum temperature alarm | AU | 9.9 | |
| 10 | Minimum temperature alarm | AL | -5.0 | |
| 11 | Temperature alarm delay | Ad | 15 | |
| 12 | Delay of temperature alarm at start up | dA | 99 | |
| 13 | THIS PARAMETER IS NOT EXPLAINED | tb | Y | |
| 14 | Evaporator probe display | d2 | nP | Read only – nP = Not Present |
| 15 | Parameter code table | Pt | 5 | Read only |
| 16 | Firmware release | rL | 1.5 | Read only |

3 Operational Matrix

| Control System Status | Dixell Status | Compressor | Unit |
|-----------------------------|--|------------------------------|---|
| Normal | Normal | Cycling under normal control | Maintaining temperature setpoint of 4C |
| Relay stuck in closed state | Normal | Cycling under Dixell control | Maintaining temperature above 2C |
| Relay stuck in open state | Normal | Off | Temperature will rise to ambient |
| Normal | Probe failure – relay opens in failsafe mode | Off | Temperature will rise to ambient |
| Normal | General power failure – relay opens on power failure | Off | Temperature will rise to ambient |
| Normal | Relay stuck open state | Off | Temperature will rise to ambient |
| Normal | Relay stuck in closed state | Cycling under normal control | Maintaining temperature |
| Relay stuck in closed state | Relay stuck in closed state | On | Temperature drops to bottom out temperature |

Dixell Operating Manual

dIXEL

Operating Manual

1598024510

DIGITAL CONTROLLER

XR02CX

1. CONTENTS

| | |
|------------------------------|---|
| 1. Contents | 1 |
| 2. General warnings | 1 |
| 3. General description | 1 |
| 4. Regulation | 1 |
| 5. Defrost | 1 |
| 6. Front panel commands | 1 |
| 7. Parameters | 1 |
| 8. Installation and mounting | 2 |
| 9. Electrical connections | 2 |
| 10. How to use the hot key | 2 |
| 11. Alarm signalling | 2 |
| 12. Technical data | 2 |
| 13. Connections | 2 |
| 14. Default setting values | 2 |

2. GENERAL WARNINGS

PLEASE READ BEFORE USING THIS MANUAL

- This manual is part of the product and should be kept near the instrument for easy and quick reference.
- The instrument shall not be used for purposes different from those described hereunder. It cannot be used as a safety device.
- Check the application limits before proceeding.

SAFETY PRECAUTIONS


- Check the supply voltage is correct before connecting the instrument.
- Do not expose to water or moisture: use the controller only within the operating limits avoiding sudden temperature changes with high atmospheric humidity to prevent formation of condensation.
- Warning: disconnect all electrical connections before any kind of maintenance.
- Fit the probe where it is not accessible by the End User. The instrument must not be opened.
- In case of failure or faulty operation send the instrument back to the distributor or to "Dixell S.p.A." (see address) with a detailed description of the fault.
- Consider the maximum current which can be applied to each relay (see Technical Data).
- Ensure that the wires for probes, loads and the power supply are separated and far enough from each other, without crossing or intertwining.
- In case of applications in industrial environments, the use of mains filters (our mod. FT1) in parallel with inductive loads could be useful.

3. GENERAL DESCRIPTION

Model XR02CX, format 32 x 74 x 50 mm, is a digital thermostat with off cycle defrost designed for refrigeration applications at normal temperature. It provides a relay output to drive the compressor. It is also provided with 2 NTC probe input. The instrument is fully configurable through special parameters that can be easily programmed through the keyboard or the by HOTKEY.

4. REGULATION

THE REGULATION OUTPUT



In case of fault in the thermostat probe the start and stop of the compressor are timed through parameters "Cy" and "Cn".

5. DEFROST

Defrost is performed through a simple stop of the compressor. Parameter "id" controls the interval between defrost cycles, while its length is controlled by parameter "Md".

6. FRONT PANEL COMMANDS



SET

To display target set point, in programming mode it selects a parameter or confirm an operation



To start a manual defrost



In programming mode it browses the parameter codes or increases the displayed value



In programming mode it browses the parameter codes or decreases the displayed value



KEYS COMBINATION

- SET + To lock or unlock the keyboard
- SET + To enter in programming mode
- SET + To return to room temperature display

| LED | MODO | SIGNIFICATO |
|-----|----------|---|
| | On | Compressore enabled |
| | Flashing | Anti short cycle delay enabled (AC parameter) |

| | | |
|--|----------|----------------------|
| | On | Defrost in progress |
| | Flashing | Dripping in progress |
| | On | Measurement unit |
| | Flashing | Programming mode |
| | On | Measurement unit |
| | Flashing | Programming mode |

HOW TO SEE THE SET POINT

- Push and immediately release the **SET** key, the set point will be showed;
- Push and immediately release the **SET** key or wait about 5s to return to normal visualisation.

HOW TO CHANGE THE SETPOINT

- Push the **SET** key for more than 2 seconds to change the Set point value;
- The value of the set point will be displayed and the "C" or "F" LED starts blinking;
- To change the Set value push the o or n arrows within 10s.
- To memorise the new set point value push the **SET** key again or wait 10s.

HOW TO START A MANUAL DEFROST (ONLY XR02CX)

Push the DEF key for more than 2 seconds and a manual defrost will start

HOW TO CHANGE A PARAMETER VALUE

To change the parameter's value operate as follows:

- Enter the Programming mode by pressing the **SET** + keys for 3s ("C" or "F" LED starts blinking).
- Select the required parameter. Press the **SET** key to display its value
- Use or to change its value.
- Press **SET** to store the new value and move to the following parameter.

To exit: Press **SET** + or wait 15s without pressing a key.

NOTE: the set value is stored even when the procedure is exited by waiting the time-out to expire.

HIDDEN MENU

The hidden menu includes all the parameters of the instrument.

HOW TO ENTER THE HIDDEN MENU

- Enter the Programming mode by pressing the **SET** + keys for 3s ("C" or "F" LED starts blinking).
- Released the keys, then push again the **SET** + keys for more than 7s. The L2 label will be displayed immediately followed from the Hy parameter.

NOW YOU ARE IN THE HIDDEN MENU.

- Select the required parameter.
- Press the **SET** key to display its value
- Use or to change its value.
- Press **SET** to store the new value and move to the following parameter.

To exit: Press **SET** + or wait 15s without pressing a key.

NOTE1: if none parameter is present in L1, after 3s the "nP" message is displayed. Keep the keys pushed till the L2 message is displayed.

NOTE2: the set value is stored even when the procedure is exited by waiting the time-out to expire.

HOW TO MOVE A PARAMETER FROM THE HIDDEN MENU TO THE FIRST LEVEL AND VICEVERSA.

Each parameter present in the HIDDEN MENU can be removed or put into "THE FIRST LEVEL" (user level) by pressing **SET** + . In HIDDEN MENU when a parameter is present in First Level the decimal point is on.

TO LOCK THE KEYBOARD

- Keep pressed for more than 3s the and keys.
- The "OF" message will be displayed and the keyboard will be locked. If a key is pressed more than 3s the "OF" message will be displayed.

TO UNLOCK THE KEYBOARD

Keep pressed together for more than 3s the and keys till the "on" message will be displayed.

7. PARAMETERS

REGULATION

Hy Differential: (0.1°C ÷ 25°C) Intervention differential for set point. Compressor Cut IN is SET POINT + differential (Hy). Compressor Cut OUT is when the temperature reaches the set point.

LS Minimum SET POINT: (-55°C+SET/-58°F+SET). Sets the minimum value for the set point.

US Maximum SET POINT: (SET+99°C/SET+99°F). Set the maximum value for set point.

ot First probe calibration: (-9.9÷9.9°C) allows to adjust possible offset of the first probe.

P2 Evaporator probe presence: n= not present; y= the defrost stops by temperature.

oE Second probe calibration: (-9.9÷9.9°C) allows to adjust possible offset of the second probe

od Outputs activation delay at start up: (0÷99min) This function is enabled at the initial start up of the instrument and inhibits any output activation for the period of time set in the parameter.

AC Anti-short cycle delay: (0÷50 min) minimum interval between the compressor stop and the following restart.

Cy Compressor ON time with faulty probe: (0÷99 min) time during which the compressor is active in case of faulty thermostat probe. With Cy=0 compressor is always OFF.

Cn Compressor OFF time with faulty probe: (0÷99 min) time during which the compressor is OFF in case of faulty thermostat probe. With Cn=0 compressor is always active.

DISPLAY

CF Measurement unit: (°C÷°F) °C =Celsius; °F =Fahrenheit. **WARNING:** When the measurement unit is changed the SET point and the values of the parameters Hy, LS, US, oE, o1, AU, AL have to be checked and modified if necessary).

rE Resolution (only for °C): (dE ÷ in) dE= decimal between -9.9 and 9.9°C; in= integer;

Ld Default display: (P1 + P2) P1= thermostat probe; P2= evaporator probe. SP=Set point

dy Display delay: (0÷15 min.) when the temperature increases, the display is updated of 1 °C/1°F after this time.

DEFROST

dE Defrost termination temperature: (-50÷50°C) if ot=Y it sets the temperature measured by the evaporator probe, which causes the end of defrost.

XR02CX US.doc

XR02CX

1/2


- id** Interval between defrost cycles: (0-99 ore) Determines the time interval between the beginning of two defrost cycles.
- md** Maximum length for defrost: (0-99 min. with 0 no defrost) when $ot=n$, (not evaporator probe: timed defrost) it sets the defrost duration, when $ot=y$ (defrost end based on temperature) it sets the maximum length for defrost.
- df** Display during defrost: (rt / it / St / dF) rt= real temperature; it= start defrost temperature; St= SET-POINT; dF= label dF.

ALARMS

- AU** Maximum temperature alarm: (AL-99°C) when this temperature is reached the alarm is enabled, after the "Ad" delay time.
- AL** Minimum temperature alarm: (-55+AU°C) when this temperature is reached the alarm is enabled, after the "Ad" delay time.
- Ad** Temperature alarm delay: (0-99 min) time interval between the detection of an alarm condition and alarm signalling.
- da** Exclusion of temperature alarm at startup: (0-99 min) time interval between the detection of the temperature alarm condition after instrument power on and alarm signalling.

OTHER

- d2** Evaporator probe display (read only)
- Pt** Parameter code table
- rL** Software release

8. INSTALLATION AND MOUNTING

Instrument **XR02CX** shall be mounted on vertical panel, in a 29x71 mm hole, and fixed using the special bracket supplied.

The temperature range allowed for correct operation is 0-60 °C. Avoid places subject to strong vibrations, corrosive gases, excessive dirt or humidity. The same recommendations apply to probes. Let air circulate by the cooling holes.

9. ELECTRICAL CONNECTIONS

The instrument is provided with screw terminal block to connect cables with a cross section up to 2.5 mm². Before connecting cables make sure the power supply complies with the instrument's requirements. Separate the probe cables from the power supply cables, from the outputs and the power connections. Do not exceed the maximum current allowed on each relay, in case of heavier loads use a suitable external relay.

9.1 PROBES

The probes shall be mounted with the bulb upwards to prevent damages due to casual liquid infiltration. It is recommended to place the thermostat probe away from air streams to correctly measure the average room temperature. Place the defrost termination probe among the evaporator fins in the coldest place, where most ice is formed, far from heaters or from the warmest place during defrost, to prevent premature defrost termination.

10. HOW TO USE THE HOT KEY**10.1 HOW TO PROGRAM THE HOT KEY FROM THE INSTRUMENT (UPLOAD)**

- Program one controller with the front keypad.
- When the controller is ON, insert the "Hot key" and push Δ key; the "uP" message appears followed by a flashing "En".
- Push "SET" key and the "En" will stop flashing.
- Turn OFF the instrument remove the "Hot Key", then turn it ON again.

NOTE: the "Er" message is displayed for failed programming. In this case push again o key if you want to restart the upload again or remove the "Hot key" to abort the operation.

10.2 HOW TO PROGRAM AN INSTRUMENT USING HOT KEY (DOWNLOAD)

- Turn OFF the instrument.
- Insert a programmed "Hot Key" into the 5 PIN receptacle and then turn the Controller ON.
- Automatically the parameter list of the "Hot Key" is downloaded into the Controller memory, the "do" message is blinking followed by a flashing "En".
- After 10 seconds the instrument will restart working with the new parameters.
- Remove the "Hot Key".

NOTE: the "Er" message is displayed for failed programming. In this case push again o key if you want to restart the upload again or remove the "Hot key" to abort the operation.

11. ALARM SIGNALLING

| Mess. | Cause | Outputs |
|-------|---------------------------|--|
| "P1" | Room probe failure | Compressor output according to "Cy" e "Cn" |
| "P2" | Evaporator probe failure | Defrost end is timed |
| "HA" | Maximum temperature alarm | Outputs unchanged |
| "LA" | Minimum temperature alarm | Outputs unchanged |
| "EA" | External alarm | Outputs unchanged |
| "CA" | Serious external alarm | All outputs OFF. |
| "dA" | Door Open | Compressor and fans restarts |

11.1 ALARM RECOVERY


Probe alarms **P1** and **P2** start some seconds after the fault in the related probe; they automatically stop some seconds after the probe restarts normal operation. Check connections before replacing the probe. Temperature alarms **HA** and **LA** automatically stop as soon as the temperature returns to normal values.

Alarms **EA** and **CA** (with IF=bL) recover as soon as the digital input is disabled.

12. TECHNICAL DATA

Housing: self extinguishing ABS.
Case: frontal 32x74 mm; depth 60mm;
Mounting: panel mounting in a 71x29mm panel cut-out
Protection: IP20; **Frontal protection:** IP65

Connections: disconnectable terminal block $\leq 2.5 \text{ mm}^2$ wiring and 6.3mm fast-on
Power supply: according to the model $\pm 10\%$; 230Vac $\pm 10\%$, 50/60Hz, 110Vac $\pm 10\%$, 50/60Hz
Power absorption: 3.5 VA max
Display: 2 digits, red LED, 14,2 mm high; **Inputs:** 2 NTC.
Relay outputs: compressor SPST 8(3) A, 250Vac; 20(8)A 250Vac
Data storing: on the non-volatile memory (EEPROM).
Kind of action: 1B; **Pollution grade:** 2; **Software class:** A;
Rated impulsive voltage: 2500V; **Overvoltage Category:** II
Operating temperature: 0-60 °C; **Storage temperature:** -30+85 °C.
Relative humidity: 20-85% (no condensing)
Measuring and regulation range: NTC -40+110°C (-40+230°F);
Resolution: 0.1 °C or 1 °C or 1 °F (selectable); **Accuracy (ambient temp. 25°C):** $\pm 0.7 \text{ °C} \pm 1 \text{ digit}$

13. CONNECTIONS

NOTE: Fast-on maximum current 16A

14. DEFAULT SETTING VALUES

| LBL | DESCRIPTION | RANGE | DEFAULT | LEVEL |
|-------------------|---|----------------------|---------|-------|
| REGULATION | | | | |
| Hy | Differential | 0.1 - 25°C/1 - 45°F | 36 °F | L1 |
| LS | Minimum Set Point | -55°C-SET/-67°F-SET | -40 °F | L2 |
| US | Maximum Set Point | SET-99°C/ SET-210°F | 99°F | L2 |
| ot | First probe calibration | -9.9-9.9°C/-18-18°F | 0.0 | L2 |
| P2 | Second probe presence | n - Y | y | L2 |
| oe | Second probe calibration | -9.9-9.9°C/-18-18°F | 0.0 | L2 |
| od | Outputs activation delay at start up | 0 + 99 min | 0 | L2 |
| AC | Anti-short cycle delay | 0 + 50 min | 0 | L1 |
| Cy | Compressor ON time faulty probe | 0 + 99 min | 15 | L2 |
| Cn | Compressor OFF time faulty probe | 0 + 99 min | 30 | L2 |
| DISPLAY | | | | |
| CF | Measurement units | °C - °F | °F | L2 |
| rE | Resolution (only for °C) | dE - in | in | L1 |
| Ld | Default Display | P1 - P2 - SP | P1 | L2 |
| dy | Display delay | 0 + 15 min | 0 | L2 |
| DEFROST | | | | |
| dE | Defrost termination temperature | -50-50°C/-58-122°F | 46 °F | L1 |
| id | Interval between defrost cycles | 0 + 99 hours | 6 | L1 |
| md | Maximum length for defrost | 0 + 99 min. | 20 | L1 |
| df | Display during defrost | rt - in - dE | it | L2 |
| ALARMS | | | | |
| AU | Maximum temperature alarm | ALL-99°C / ALL-210°F | 99 °F | L2 |
| AL | Minimum temperature alarm | -55°C-ALU/-67°F-ALU | -50 °F | L2 |
| Ad | Temperature alarm delay | 0 + 99 min | 15 | L2 |
| da | Exclusion of temperature alarm at startup | 0 + 99 min | 99 | L2 |
| OTHER | | | | |
| d2 | Evaporator probe display | Read Only | --- | L1 |
| Pt | Parameter code table | Read Only | --- | L2 |
| rL | Firmware release | Read Only | --- | L2 |

dIXEL S.p.a.


Z.I. Via dell'Industria, 27 - 32010 Pieve d'Alpago (BL) ITALY
 tel. +39 - 0437 - 98 33 - fax +39 - 0437 - 98 93 13
<http://www.dixel.com> E-mail: dixel@dixel.com

Anti-Tip Bracket Instructions



WARNING: SAFETY INSTRUCTIONS

- REFRIGERATOR/FREEZER MUST BE SECURED BY THE ANTI-TIP BRACKET SUPPLIED IN THIS BAG.
- UNLESS PROPERLY INSTALLED, REFRIGERATOR/FREEZER COULD TIP WHEN SHELVES/ DRAWERS ARE LOADED. INJURY AND DAMAGE TO EQUIPMENT AND CONTENTS MAY RESULT FROM REFRIGERATOR/ FREEZER TIPPING
- THIS REFRIGERATOR/FREEZER HAS BEEN DESIGNED TO MEET ALL RECOGNIZED INDUSTRY TIP STANDARDS FOR ALL NORMAL CONDITIONS.




INSTALLATION INSTRUCTIONS:

Installation instructions are provided for wood and concrete floors. Any other type of construction may require special installation techniques as deemed necessary to provide adequate fastening of the Anti-Tip bracket to the floor. For installation on floors other than wood or concrete, please contact technical support.

The use of this bracket does not prevent the tipping of the Refrigerator/Freezer when not properly installed.

Materials Supplied

1. Anti-Tip Bracket (1)
2. 5/16" Lag Bolt (2)
3. Lag Screw Anchor (2), for concrete installation only
4. Bracket location template



Tools Required

Wood Floor

- Flashlight
- Tape Measure
- Drill
- 15/64" (6mm) Drill Bit
- 1/2" (13mm) Wrench
- 3/4" (19mm) Wrench

Concrete Floor

- Flashlight
- Tape Measure
- Hammer Drill
- 1/2" (13mm) Masonry Bit
- 1/2" (13mm) Wrench
- 3/4" (19mm) Wrench

Step 1

Locating the Bracket

- a. Determine where you want the centerline of the refrigerator/freezer to be
- b. Place the included template on the floor lined up with the centerline of the refrigerator/freezer and keep 6"-12" between the wall and the back of the unit
- c. On the floor, mark the location of Hole #1 & Hole #2 (also Hole #3 & Hole #4 for 50ft³ & 75ft³ models).

Step 2

Anti-Tip Bracket Installation

Wood Construction

- a. Drill 15/64" (6mm) pilot holes in locations marked in step 1
- b. Place bracket on floor aligned with holes
- c. Use supplied lag bolts to attach bracket to floor

Concrete Construction

- a. Drill 1/2" (13mm) holes in locations marked in step 1 with masonry bit
- b. Slide Lag Screw Anchors into holes to be flush with floor surface
- c. Place bracket on floor aligned with holes
- d. Use supplied lag bolts to attach bracket to floor

Step 3

Adjusting Bolt in Refrigerator/Freezer

- a. Locate 1/2" bolt attached to bottom of cabinet
- b. Unscrew 1/2" bolt until there is the required clearance between floor and head of bolt as shown in Figure 1
- c. Tighten lock nut against bottom of unit

Step 4

Refrigerator/Freezer Positioning

- a. Line up 1/2" bolt installed in Step 3 with anti-tip bracket
- b. Roll or slide Refrigerator/Freezer into position until bolt stops against bracket
- c. Lock the casters

Step 5

Checking the Installation

- a. Complete the installation of the Refrigerator/Freezer per the installation instructions provided with the product.
- b. Check to see if the Anti-Tip bracket is installed properly by shining light under cabinet and confirming bolt in cabinet is secured by bracket on floor




FIGURE 1

4-20mA Temperature Transmitter Information

Testing of Self Powered 4-20mA Temperature Transmitter Options Option #'s 7203-AV, 7203FMS, 7204, 7205TA (or 4-20mA loops in general)

Main parts – 326518H01 through H04 TempTrans (Minco TT518), 327624H01 24v DC Power Supply or 16.5v DC Power Supply and a dummy load resistor (1/4w 5%) – 25 ohms or 250 ohms
Refer to instructional drawing #s 311977I02, 325436H01 FMS and wiring diagram #s 325439H01 FMS, 326861H01 General, and 327987H01 TSX.

Theory in general –

The power supply provides 24v DC to the loop. The acceptable voltage range is 8 to 30v DC for customers who supply their own power to the loop. The Minco TempTrans converts the sensor (RTD) reading into a milli Amp signal based on Ohm's Law ($V=IR$). 4mA represents the coldest end of the scale printed on the back of the TempTrans, and 20mA represents the warmest end of the scale.

| | | | |
|-----|--------------------------------|--------------|---------------|
| For | 326518H01 (Frzs -20/-30/-86s), | 4mA = -100°C | & 20mA = 0°C |
| | 326518H02 (Refrigerators), | 4mA = 0°C | & 20mA = 40°C |
| | 326518H04 (FMS 2320, & 2305), | 4mA = -66°C | & 20mA = 54°C |

Testing -

Testing can be done by current or voltage measurements. A dummy load resistor will need to be added temporarily across the terminal strip connections "TempTrans Neg" and "24v Supply Neg" or "16.5 Supply Neg" (see 311977I02 or 327987H01). The resistor simulates the customers test equipment (see wiring diagrams for "control room - R1"). The size of the resistor is limited based on maximum resistance of the loop formula:

$$R \text{ loop max} = (V \text{ supply} - 8v) / 0.023 \text{ amps, in our 24v case, } R \text{ loop max} = 696 \text{ ohms}$$

Note: for HPLRF TSX, the supply is 16.5v, so $R \text{ loop max} = 370 \text{ ohms}$

Measurements only need to reflect that the transmitter is functioning near the expected range since actual readings will be dependent on additional factors at the customer's site. The best time to test is while the cabinet is operating within its TempTrans range, but for cabinets that are out of range (-86 at ambient), the TempTrans will show the "out of range" value of 3.5mA if the circuit is connected correctly.

Current Method

For current measurements, the size of the resistor just needs to be less than R loop max or the resistor can be replaced at the terminal strip by a DC ammeter. Break into the circuit with an Amp meter set to the mA DC scale. The current reading from the meter can be compared to the calculated expected value to confirm that it corresponds to the cabinet temperature. To calculate the expected value, use the formula:

Current for T_{sensor} = $((T_{\text{sensor}} - T_{\text{zero}} / T_{\text{span}} - T_{\text{zero}}) * 16\text{mA}) + 4\text{mA}$ T_{sensor} = the temp at the sensor

T_{zero} = the printed temp for 4mA (H02 = 0°C)

T_{span} = the printed temp for 20mA (H02 = 40°)

For example: If a cabinet using 38197H02 is reading +4°C Then 4°C =

$$((4-0 / 40-0)*16) + 4 = 5.6\text{mA}$$

Or

A cabinet using H01 is reading -86°C

$$\text{Then } -86^{\circ}\text{C} = ((-86+100 / 0+100)*16) + 4 = 6.24\text{mA}$$

Reduced °C formulas: for H01, current at T_{sensor} = $(T_{\text{sensor}} * .16) + 4\text{mA}$
 for H02, current at T_{sensor} = $(T_{\text{sensor}} * .4) + 4\text{mA}$
 for H04, current at T_{sensor} = $(T_{\text{sensor}} * .1333) + 4\text{mA}$

Voltage Method

For voltage measurements, the size of the resistor effects the voltage range. For example: a 250 ohm resistor will cause a voltage range of 1v to 5v DC, and a 25 ohm resistor will cause a range of .1v to .5v DC. Other values up to R loop max are acceptable but require recalculating with Ohms Law ($V=IR$) to generate the expected voltage range.

Place the meter leads across the dummy load resistor to read the voltage. To calculate the expected voltage, use the calculated current value times the resistor value.

For a 250 ohm resistor with H02, the voltage at the dummy load resistor of a refrigerator at

$$+4^{\circ}\text{C} = 5.6\text{mA} * 250 = 1.4\text{v}$$

and using a 250 ohm resistor with H01, the voltage at the dummy load resistor of a ULT at

$$-86^{\circ}\text{C} = 6.24\text{mA} * 250 = 1.56\text{v}$$

You can also use percentage approximation to tell if the TempTrans is working. Based on a 250 ohm resistor and the 1v to 5v range, temps at 50% of the range will read ~ 3v, 25% will read ~2v, 75% will read ~ 4v, etc. For example, a refrigerator at ambient will be in the 3.0v to 3.5v area, since 20°C = 3v and 25°C = 3.5v; and a refrigerator at normal operating temp should be in the 1.0v to 1.5v range since 0°C = 1v and 10°C = 1.5v.

Quick Reference Charts are available mapping out the TempTrans ranges with both 25 and 250 ohm resistors as are Resistance vs. Temperature charts for 100 PRTD with .385 coefficient.


Possible problems – other than a bad component, open/shorted circuits (low voltage or line voltage), reversed polarity on the power supply circuit, or R loop max being exceeded by the customer -

- 1) Power Supply accepts 100v to 240v AC as an input. Output should be 24v DC and the green LED should be on. In the 16.5v case, the TempTrans is piggybacking off of the units supply, so there is no green LED
- 2) The voltage across the TempTrans should read 24v (16.5) minus the voltage drop of the dummy load resistor (building system). When the TempTrans sees less than 8v, it will not be able to regulate the loop at the 20mA end, so mA values will start dropping off.
- 3) Reversing the 24v wires at the power supply or the TempTrans causes the voltage across the dummy load to read 0v.

Additional resource – www.minco.com

4-20mA Option Numbers

Non-RoHS Temp Transmitter Options (Blue)




| 4-20mA Thermo P/N (Transmitter Only) | Temp Range @ 4mA | Temp Range @ 20mA | Application | Loop Powered Transmitter and Sensor Option Number | Self Powered Transmitter and Sensor Option Number | Field install Loop Powerd | Field install Self Powerd | Minco Model Number |
|---|---------------------|----------------------|--------------------------------|--|---|------------------------------|------------------------------|--------------------|
| 38197H01 | -100°C | 0.0°C | ULT Freezer | 4706 | 7203-AV | | | TT111PD1EZ |
| 38197H01 | -40°C | 0.0°C | -20°C / -30°C HPLRF Freezer | 4706 | 7205TA | | | TT111PD1EZ |
| 38197H02 | 0.0°C | 40°C | +4°C Refrigerator | 6907-1 | 7204 | | | TT111PD1EG |
| 38197H02 | 0.0°C | 40°C | +4°C Refrigerator | 6907-2 | 7204 | 6907-2 | | TT111PD1EG |
| 38197H03 | -185°C | -95°C | -140 / -150 Freezers | 7261 | N/A | | | TT111PD1KE |
| 38197H04** | -66°C | 54°C | FMS Refrigerator / Freezers | 4706FMS | 7203FMS | | | TT111PD1DB |
| 38197H04 | -20°C | 100°C | BOD's / Incubators | 7258 | 7259 | | | TT111PD1DB |


RoHS Temp Transmitter Options (Red Puck) July 2014 to Current




| 4-20mA Thermo P/N (Transmitter Only) | Temp Range @ 4mA | Temp Range @ 20mA | Application | Loop Powered Transmitter and Sensor Option Number | Self Powered Transmitter and Sensor Option Number | Field install Loop Powerd | Field install Self Powerd | Minco Model Number |
|---|---------------------|----------------------|--------------------------------|--|---|------------------------------|------------------------------|---------------------|
| 326518H01 | -100°C | 0.0°C | ULT Freezer | 4706 | 7203-AV | 7209TA | 7207TA | TT518PD(-100/0)C1Y |
| 326518H01 | -40°C | 0.0°C | -20°C / -30°C HPLRF Freezer | 4706 | 7205TA | | | TT518PD(-100/0)C1Y |
| 326518H02 | 0.0°C | 40°C | +4°C Refrigerator | 6907-1 | 7204 | 7208TA | 7206TA | TT518PD(0/40)C1Y |
| 326518H02 | 0.0°C | 40°C | +4°C Refrigerator | 6907-2 | | | | TT518PD(0/40)C1Y |
| 326518H04** | -66°C | 54°C | FMS Refrigerator / Freezers | 4706FMS | 7203FMS | | | TT518PD(-20/100)C1Y |

** Offset due to use of intrinsically safe barrier with H04 Transmitter







NOTES:

1. FOR WIRE HARNESS CONSTRUCTION & SPECIFICATIONS SEE TFS 28305H01.
2. COMPONENTS AND WIRE ARE AS SPECIFIED OR AN EQUIVALENT APPROVED BY THERMOFISHER ENGINEERING.
3. OPTIONS SHOWN WITH DASHED LINES.
4. ALL DARK LINES CONTAIN HIGH POTENTIAL VOLTAGE.

ARTICLE OR MATERIAL MUST COMPLY WITH THE ENVIRONMENTAL COMPLIANCE SPECIFICATION TFS DRAWING 260111S01

| | | | |
|-----|--------|------------|--|
| E | 005332 | 2019-04-04 | UPDATED 12.30 CU FT HARNESS DETAILS AND PEP HEATER HARNESS |
| D | 005511 | 2019-02-25 | ADDED FRENCH NOTES IN SHEET-2 |
| C | 005144 | 2018-10-11 | ADDED SHEET 2 WITH THE FRENCH LANGUAGE |
| B | 004960 | 2018-07-24 | UPDATED POWER SWITCH AND POWER INLET HARNESS |
| A | 004610 | 2018-03-06 | INITIAL RELEASE |
| REV | DCO # | DATE | DESCRIPTION |

| | | | |
|--|--------------|--|---|
| REVISIONS | | | THE INFORMATION CONTAINED HEREIN IS THE PROPRIETARY DATA OF THERMO FISHER SCIENTIFIC. IT IS NOT TO BE REPRODUCED OR DISCLOSED IN WHOLE OR IN PART WITHOUT PRIOR WRITTEN CONSENT OF THERMO FISHER SCIENTIFIC. DO NOT SCALE THIS DRAWING. |
| ThermoFisher S C I E N T I F I C | | | |
| 275 Aiken Road Asheville, NC 28804 | | | |
| TITLE DIAGRAM, ELECTRICAL WIRING -20 CW | | | |
| 3rd ANGLE PROJECTION | | | DWG. NO. 328059H06 |
| SIZE C | SHEET 1 OF 2 | | |

Charge Information

| TSX Charge Info | | | | | |
|-----------------|-------------------|-----------|-------|-------|------|
| Size CuFt | Suffix | Temp °C | Lbs | Grams | Oz |
| 12 | 17 | +4° / +5° | 0.242 | 110 | 3.88 |
| 12 | 17 | -30° | 0.209 | 95 | 3.35 |
| 23 | 16 / 17 | +4° / +5° | 0.319 | 145 | 5.11 |
| 23 | 18 | +4° / +5° | 0.308 | 140 | 4.94 |
| 23 | 16 / 17 / 18 / 19 | -20° | 0.275 | 125 | 4.41 |
| 23 | 20 | -20° | 0.319 | 145 | 5.11 |
| 23 | 16 | -30° | 0.253 | 115 | 4.06 |
| 23 | 17 / 18 | -30° | 0.242 | 110 | 3.88 |
| 30 | 17 | +4° / +5° | 0.308 | 140 | 4.94 |
| 30 | 18 / 19 | -20° | 0.319 | 145 | 5.11 |
| 30 | 17, 18 | -30° | 0.286 | 130 | 4.59 |
| 45 | 17 | +4° / +5° | 0.319 | 145 | 5.11 |
| 50 | 17 | +4° / +5° | 0.319 | 145 | 5.11 |
| 50 | 17 | -30° | 0.297 | 135 | 4.76 |

Pressure and Amp Draw

TSX 23FT Model Pressure & Current Data (20C ambient cycle)

| Model | Suffix | P1max(psig) | P1min(psig) | P2max(psig) | P2min(psig) | Unit Amps (compressor on) | Unit Amps (compressor off) |
|--------------------|---------|-------------|-------------|-------------|-------------|---------------------------|----------------------------|
| TSX2305A, TSX2304A | 16 - 17 | 160 | 155 | 70.0 | 27.0 | 1.84 | 0.97 |
| TSX2305D, TSX2304D | 16 - 17 | 160 | 155 | 70.0 | 27.0 | 1.60 | 0.68 |
| TSX2305V, TSX2304V | 16 - 17 | 160 | 155 | 70.0 | 27.0 | 1.60 | 0.80 |
| TSX2305A, TSX2304A | 18 | | | | | | |
| TSX2305D, TSX2304D | 18 | | | | | | |
| TSX2305V, TSX2304V | 18 | | | | | | |
| TSX2320A | 16-19 | 145 | 130 | 15.0 | 5.3 | 3.00 | 1.30 |
| TSX2320D | 16-19 | 145 | 130 | 15.0 | 5.3 | 1.60 | 0.55 |
| TSX2320V | 16-19 | 145 | 130 | 15.0 | 5.3 | 1.50 | 0.58 |
| TSX2320A | 20 | 125 | 120 | 6.7 | 3.7 | 3.00 | 1.30 |
| TSX2320D | 20 | 125 | 120 | 6.7 | 3.7 | 1.60 | 0.55 |
| TSX2320V | 20 | 125 | 120 | 6.7 | 3.7 | 1.50 | 0.58 |
| TSX2330A | 16 | | | | | | |
| TSX2330D | 16 | | | | | | |
| TSX2330V | 16 | | | | | | |
| TSX2330A | 17, 18 | 140 | 137 | 7.0 | 2.5 | 4.50 | 1.40 |
| TSX2330D | 17, 18 | 140 | 137 | 7.0 | 2.5 | 2.20 | 0.57 |
| TSX2330V | 17, 18 | 140 | 137 | 7.0 | 2.5 | 2.15 | 0.60 |

TSX 45~50FT Model Pressure & Current Data (20C ambient cycle)

| Model | Suffix | P1max(psig) | P1min(psig) | P2max(psig) | P2min(psig) | Unit Amps(compressor on) | Unit Amps(compressor off) |
|--------------------|--------|-------------|-------------|-------------|-------------|--------------------------|---------------------------|
| TSX5030A | 17 | 145 | 138 | 10 | 4 | 8.70 | 1.80 |
| TSX5030D | 17 | 145 | 138 | 10 | 4 | 4.80 | 0.80 |
| TSX5030V | 17 | 145 | 138 | 10 | 4 | 4.40 | 1.00 |
| TSX4505A | 17 | 144 | 142 | 67 | 34 | 3.50 | 1.60 |
| TSX4505D | 17 | 144 | 142 | 67 | 34 | 1.90 | 1.00 |
| TSX4505V | 17 | 144 | 142 | 67 | 34 | 1.40 | 0.80 |
| TSX5005A, TSX5004A | 17 | 144 | 142 | 67 | 34 | 4.10 | 2.50 |
| TSX5005D, TSX5004D | 17 | 144 | 142 | 67 | 34 | 1.70 | 0.95 |
| TSX5005V, TSX5004V | 17 | 144 | 142 | 67 | 34 | 1.60 | 0.80 |

TSX 12FT Model Pressure & Current Data (20C ambient cycle)


| Model | Suffix | P1max(psig) | P1min(psig) | P2max(psig) | P2min(psig) | Unit Amps(compressor on) | Unit Amps(compressor off) |
|--------------------|--------|-------------|-------------|-------------|-------------|--------------------------|---------------------------|
| TSX1205A, TSX1204A | 17 | 152 | 143 | 64 | 38 | 2.30 | 1.40 |
| TSX1205D, TSX1204D | 17 | 152 | 143 | 64 | 38 | 1.40 | 1.00 |
| TSX1205V, TSX1204V | 17 | 152 | 143 | 64 | 38 | 1.30 | 0.95 |
| TSX1230A | 17 | 137 | 134 | 6.8 | 2.5 | 3.30 | 1.00 |
| TSX1230D | 17 | 137 | 134 | 6.8 | 2.5 | 1.65 | 0.60 |
| TSX1230V | 17 | 137 | 134 | 6.8 | 2.5 | 1.70 | 0.60 |

TSX 30FT Model Pressure & Current Data (20C ambient cycle)

| Model | Suffix | P1max(psig) | P1min(psig) | P2max(psig) | P2min(psig) | Unit Amps(compressor on) | Unit Amps(compressor off) |
|--------------------|---------|-------------|-------------|-------------|-------------|--------------------------|---------------------------|
| TSX3005A, TSX3004A | 17 | 134 | 130 | 65.3 | 37.3 | 3.50 | 1.70 |
| TSX3005D, TSX3004D | 17 | 134 | 130 | 65.3 | 37.3 | 1.60 | 0.80 |
| TSX3005V, TSX3004V | 17 | 134 | 130 | 65.3 | 37.3 | 1.60 | 0.90 |
| TSX3020A | 18 - 19 | 128 | 123 | 13.3 | 4.8 | 2.70 | 1.10 |
| TSX3020D | 18 - 19 | 128 | 123 | 13.3 | 4.8 | 1.60 | 0.60 |
| TSX3020V | 18 - 19 | 128 | 123 | 13.3 | 4.8 | 1.50 | 0.60 |
| TSX3030A | 17 | 135 | 126 | 6.8 | 2.5 | 3.90 | 1.20 |
| TSX3030D | 17 | 135 | 126 | 6.8 | 2.5 | 1.70 | 0.60 |
| TSX3030V | 17 | 135 | 126 | 6.8 | 2.5 | 2.00 | 0.70 |
| TSX3030A | 18 | | | | | | |
| TSX3030D | 18 | | | | | | |
| TSX3030V | 18 | | | | | | |


Notes: P1max/P2max are maximum discharge/suction pressure when compressor is on(not in defrost cycle).
P1min/P2min are minimum discharge/suction pressure when compressor is on(not in defrost cycle).
Unit Amps data are collected when door frame heaters are on (except for 4505 models which are without door frame heaters)
Data are estimated pressure/current and for reference only.

| S. NO. | CATEGORY | 12 Cu ft REF P/N | | | 23 Cu ft REF P/N | | | 30 Cu ft REF P/N | | | 45/50 Cu ft REF P/N | | |
|--------|---------------------------|------------------|-----------|---------------|------------------|-----------|---------------|------------------|-----------|---------------|---------------------|-----------|---------------|
| | | 120V/60Hz | 230V/50Hz | 208-230V/60Hz | 120V/60Hz | 230V/50Hz | 208-230V/60Hz | 120V/60Hz | 230V/50Hz | 208-230V/60Hz | 120V/60Hz | 230V/50Hz | 208-230V/60Hz |
| 1 | COMPRESSOR ASSY | 327001G05 | | | 327001G04 | | | | | | | | |
| 2 | INVERTER ASSY | 327002H07 | 327002H08 | | 327100G08 | 327100G09 | | 327100G08 | 327100G09 | | 327100G08 | 327100G09 | |
| 3 | CONDENSER | 328321H02 | | | 328321H01 | | | | | | | | |
| 4 | CONDENSER FAN ASSY | 327005G02 | 327005G01 | | 327005G02 | 327005G01 | | 327005G02 | 327005G01 | | 327005G02 | 327005G01 | |
| 5 | FILTER/DRIER | 75380H01 | | | | | | | | | | | |
| 6 | RESTRICTOR | 47417H106 | | | | | | | | | | | |
| 7 | EVAPORATOR ASSY | 29333G28 | 29333G29 | | 29333G22 | 29333G23 | | 29333G26 | 29333G27 | | 60595H33 | 60595H34 | |
| | ALTERNATE EVAPORATOR ASSY | NA | | | | | | | | | 60595H43 | 60595H44 | |
| 8 | EVAPORATOR FAN | 313261H01 | | | | | | | | | 306637H01 | 307462H03 | |
| | ALTERNATE EVAPORATOR FAN | NA | | | | | | | | | 39693H01 | 39693H02 | |
| 9 | PERIMETER HEATER | 329588H01 | 329588H02 | | 329588H03 | 329588H04 | | 329588H05 | 329588H06 | | 329588H03 | 329588H04 | |
| 10 | CONTROL PROBE | 326737H03 | | | | | | | | | | | |
| 11 | DISPLAY PROBE | 309159H09 | | | | | | | | | 309159H10 | | |
| 12 | EVAP COIL DEFROST SENSOR | 326737H05 | | | | | | | | | | | |
| 13 | AMBIENT PROBE | 326739H01 | | | | | | | | | | | |
| 14 | CONDENSATE PAN HEATER | 328840H01 | 328840H02 | | 328840H01 | 328840H02 | | 328840H01 | 328840H02 | | 328840H01 | 328840H02 | |



REF, 45/50 Cu ft: SHOWN
12, 23, 30 Cu ft: FOIL TAPE HEAT EXCHANGE
REPLACES 2" SOLDER HEAT EXCHANGE
*NO CONDENSATE PAN LOOP AND HEATER
FOR 12 Cu ft & 30 Cu ft

| S. NO. | CATEGORY | 12 Cu ft FRZ P/N | | | 23 Cu ft FRZ P/N | | | 30 Cu ft FRZ P/N | | | 50 Cu ft FRZ P/N | | |
|--------|------------------------------|------------------|-----------|---------------|------------------|-----------|----------------|------------------|-----------|---------------|------------------|-----------|---------------|
| | | 120V/60Hz | 230V/50Hz | 208-230V/60Hz | 120V/60Hz | 230V/50Hz | 208- 230V/60Hz | 120V/60Hz | 230V/50Hz | 208-230V/60Hz | 120V/60Hz | 230V/50Hz | 208-230V/60Hz |
| 1 | COMPRESSOR ASSY | 327001G04 | | | 327001H06 | | | 327001G03 | | | | | |
| 2 | INVERTER | 327100G08 | 327100G09 | | 327002H09 | | | 327100G06 | 327100G07 | | 327100G06 | 327100G07 | |
| 3 | CONDENSER | 328321H02 | | | GT350222H01 | | | 328321H01 | | | | | |
| 4 | CONDENSER FAN ASSY | 327005G02 | 327005G01 | | 327005G02 | 327005G01 | | 327005G02 | 327005G01 | | 327005G02 | 327005G01 | |
| 5 | FILTER/DRIER | 75380H01 | | | | | | | | | | | |
| 6 | RESTRICTOR | 47417H107 | | | | | | | | | 87812H07 | | |
| 7 | EVAPORATOR ASSY | 60595H37 | 60595H38 | | 327738G03 | 327738G04 | | 327738G03 | 327738G04 | | 60595H35 | 60595H36 | |
| | ALTERNATE EVAPORATOR ASSY | NA | NA | | NA | NA | | 60595H39 | 60595H40 | | 60595H41 | 60595H42 | |
| 8 | EVAPORATOR FAN | 39693H01 | 39693H02 | | 306637H01 | 307462H03 | | 306637H01 | 307462H03 | | 306637H01 | 307462H03 | |
| | ALTERNATE EVAPORATOR FAN | NA | NA | | NA | NA | | 39693H01 | 39693H02 | | 39693H01 | 39693H02 | |
| 9 | DRAIN PAN HEATER | 328269H05 | 328269H06 | | 328269H03 | 328269H04 | | 328269H03 | 328269H04 | | 328269H01 | 328269H02 | |
| 10 | DRAIN TUBE HEATER | 311730H01 | 311730H02 | | 311730H01 | 311730H02 | | 311730H01 | 311730H02 | | 311730H01 | 311730H02 | |
| 11 | PERIMETER HEATER | 329588H01 | 329588H02 | | 329588H03 | 329588H04 | | 329588H05 | 329588H06 | | 329588H03 | 329588H04 | |
| 12 | CONTROL PROBE | 326737H03 | | | | | | | | | | | |
| 13 | DISPLAY PROBE | 309159H09 | | | | | | | | | 309159H10 | | |
| 14 | EVAP COIL DEFROST SENSOR | 326737H05 | | | | | | | | | | | |
| 15 | DRAIN PAN DEFROST SENSOR | 326739H03 | | | | | | | | | | | |
| 16 | AMBIENT PROBE | 326739H01 | | | | | | | | | | | |
| 17 | SAFETY THERMOSTAT | 327883H01 | | | | | | | | | | | |
| 18 | ACCUMULATOR | 300367H01 | | | | | | | | | | | |
| 19 | HOT GAS DEFROST VALVE | 328498H01 | | | | | | | | | | | |
| 20 | HOT GAS DEFROST SOLENOID | 328499H01 | 328499H02 | | 328499H01 | 328499H02 | | 328499H01 | 328499H02 | | 328499H01 | 328499H02 | |
| 21 | TRANSFORMER | NA | | | GT350729H01 | NA | | | | | | | |
| 22 | HIGH PRESSURE CUT OFF SWITCH | NA | | | | | | 312598H04 | | | | | |




FRZ, 50 Cu ft: SHOWN
 12, 23, 30 Cu ft: FOIL TAPE HEAT EXCHANGE
 REPLACES 2" SOLDER HEAT EXCHANGE

Touch Up Paint Information

- New TSX look with new colors

- 45200H105 Facet Light Grey, Aerosol
- 45200H106 Facet Light Grey, Liquid
- 45200H108 Facet Grey, Aerosol
- 45200H109 Facet Grey, Liquid



Technical Service Bulletins of Note

| Description |
|---|
| 1708-REF-AVL-0500 TSX Distribution board connections |
| 1807-REF-AVL-0538 Launch of the TSX 45 & 50 ft units |
| 1905-REF-AVL-0569 Launch of the TSX 12, 30 ft units |
| 1909-REF-AVL-0579 -30C Defrost Disabled issue |
| 2001-REF-AVL-0589 5030 defrost icing issues with note |
| 2004-REF-AVL-0600 Audit Phase software release and new -20 config |
| 2007-REF-AVL-0607 HPLRF 12ft missing counterweights |
| 2011-REF-AVL-0616 TSX -30 defrost icing issues – 12, 23, & 30' |
| 2106-REF-AVL-0634 2320HA FMS launch, and 2320 standardization |
| 2203-REF-AVL-0669 Cardinal Anthony to CIG Glass Door Replacement |
| 2202-REF-AVL-0664 TSX Anthony to CIG Glass Door Replacement |
| 2205-REF-AVL-0675 TSX2304-05 and CH2304-05 Standardization Suffix 18 |
| 2206-REF-AVL-0681 TSX3030 Suffix 18 Change Heatcraft Evaporator to Peerless |
| 2207-REF-AVL-0687 TSX5030 Shock Loop Compressor Version 2 |

| Description |
|---|
| 326300G94 Kit instructions for -30 icing issues |
| 326300G96 Kit instructions for normal replacement board kit |
| TSX Display offset instructions |
| TSX Control offset instructions |
| HPLRF Master Parts list |
| HPLRF Parameter Tables with notes and examples |

Available Service Kits

| Compressor / Inverter Service Kits | |
|------------------------------------|---|
| Kit Number | Description |
| 303255G47 | Service Assembly, Compressor, (VNEK207U, oil & filter dryer) 23' ref |
| 303255G49 | Service Assembly, Compressor, (VNEU217U, oil & filter dryer) 30/50' freezer |
| 303255G50 | Service Assembly, Compressor, (VEGT8U & filter dryer) 23/30/45/50' ref & 1230 and 3020 freezers |
| 303255G51 | Service Assembly, Compressor, (FMFT406U, oil & filter dryer) 12' ref |
| 303255G52 | Service Assembly, Compressor, (DSL - VNEU217U VER 2) |

| 23' -30 & -20 HP to Maia Compressor Kits | | |
|--|---|--|
| Kit Number | Use | Kit Description |
| 100793G01 | TSX2330 Suffix 16 A/D | TSX2330 Inver/Comp Suffix 16 A/D SER UPG Kit |
| 100794G01 | TSX2330 Non-Suffix 16 A w/327002H01 Inverter & 327001H01 Compressor | TSX2330 Inver/Comp A SER UPG Kit |
| 100795G01 | TSX2330 Non-Suffix 16 D w/327002H01 Inverter & 327001H01 Compressor | TSX2330 Inver/Comp D/V SER UPG Kit |
| 327727G11 S | TSX2330 Non-Suffix 16 A w/327002H09 Inverter & 327001H06 Compressor | SRV ASY, COND UNIT FRZ A -30 HG |
| 327727G12 S | TSX2330 Non-Suffix 16 D w/327002H09 Inverter & 327001H06 Compressor | SRV ASY, COND UNIT FRZ D/V -30 |
| 100796G01 | TSX2320 Non-Suffix 20 D/V Volt | TSX2320 Inver/Comp A SER UPG Kit |
| 100797G01 | TSX2320 Non-Suffix 20 A Volt | TSX2320 Inver/Comp D/V SER UPG Kit |
| 327727G07 S | TSX2320 Suffix 20 A | ASSY, COND UNIT, REF A, 45/50 |
| 327727G08 S | TSX2320 Suffix 20 D/V | SRV ASY, COND UNIT REF D/V 45/ |
| 100798G01 | TSX2330 Non-Suffix 16 A w/327002H09 Inverter & 327001H06 Compressor | SRV ASY, COMPR, VNEK213U MAIA |

| Inverter Service Kits | | |
|-----------------------|-----------------------------------|-----------------------------|
| Kit Number | Description | Use |
| 327100G04 | ASSEMBLY, INVERTER, VNEK207U 115V | 23' Refrig |
| 327100G05 | ASSEMBLY, INVERTER, VNEK207U 230V | 23' Refrig |
| 327100G06 S | ASSEMBLY, INVERTER, VNEU217U 115V | 30 / 50' Frz |
| 327100G07 S | ASSEMBLY, INVERTER, VNEU217U 230V | 30 / 50' Frz |
| 327100G08 S | ASSEMBLY, INVERTER, VEGT8U 115V | 30/45/50' Refrig, 1230 3020 |
| 327100G09 S | ASSEMBLY, INVERTER, VEGT8U 230V | 30/45/50' Refrig, 1230 3020 |

| 23' -30 & -20 HP to Maia Inverter Kits | | |
|--|---|----------------------------------|
| Kit Number | Use | Kit Description |
| 100799G01 | TSX2330 Non-Suffix 16 A w/327002H09 Inverter & 327001H06 Compressor | SRV INVERTER ASY, MAIA 213u 115V |
| 100800G01 | TSX2330 Non-Suffix 16 A w/327002H09 Inverter & 327001H06 Compressor | SRV INVERTER ASY, MAIA 213u 230V |

| Heatcraft to Peerless Upgrade Kits | | |
|------------------------------------|-------------------------------|--|
| Kit Number | Use | Kit Description |
| 100807G01 | UPGD HC to Peerless EVAP 115v | Upgrade from Heatcraft to Peerless Evap. TSX3030 |
| 100808G01 | UPGD HC to Peerless EVAP D/V | Upgrade from Heatcraft to Peerless Evap. TSX3030 |

| Replacement CPU kits with the appropriate Gxx EPROM and matching Relay board code – | |
|---|---|
| Kit Number | Description |
| 326300G94 | PCB SERV KIT, TSX CPU/G57 and H02 Relay Board |
| 326300G96 | PCB SERV KIT, TSX CPU/G54 and H02 Relay Board |

For a complete list of Service Parts, please reference the Knowledge Base for related articles including #330015 Parts List.

Thermo Fisher Scientific
275 Aiken Road
Asheville, NC 28804 U.S.A.
866.984.3766
ALL RIGHTS RESERVED