


Thermo Scientific Midi 40 Small Capacity CO₂ Incubator



Technical Specifications		
Dimensions		
Exterior (w x h x d) inch/cm	18.5 x 18 x 23.5 (47 x 46.5 x 59.7)	
Interior (w x h x d) inch/cm	12 x 14 x 14 (30.5 x 35.5 x 35.5)	
Construction		
Cabinet configuration	bench top	
Interior	non-corrosive stainless steel, type 304, mirror finish	
Interior volume	1.4 cubic feet (39.6 L)	
Exterior	18 gauge, cold rolled steel, high impact painted finish	
Shelves		
Construction	stainless steel, perforated	
Shelf dimensions (w x d) inch/cm	13.5" x 11.5" (34.3 x 29.2)	
Surface area per shelf	1.08 sq. ft. (0.1 sq m)	
Temperature		
Range	5°C above ambient to 60°C	
Control accuracy	± 0.1°C	
Uniformity	± 0.5°C @ 37°C	
Alarm	user programmable high/low	
CO ₂		
Range	0-20%	
Control accuracy	± 0.1%	
Sensor	long life thermal conductivity sensor	
Uniformity	0.10%	
Inlet pressure required	15 PSIG (1.0 bar), 99.5 medical grade CO ₂	
Inlet air filtration	0.3 micron antimicrobial filter	
Alarm	user programmable high/low	
Humidity		
rH	to 95% @ 37°C	
Source	humidity pan, removable	
Electrical		
Rated voltage	120V, 60 Hz	230V, 50/60 Hz
Amps/breaker	1.5/15	1.0/15
Weight		
Net	60 lbs (27 kg)	
Shipping	70 lbs (32 kg)	
Warranty	1 year	

All specifications subject to change without notice.

Ordering Information	
Cat. No.	Description
3403	Thermo Scientific Midi 40 Small Capacity CO ₂ Incubator, 120 V, 60 Hz
3404	Thermo Scientific Midi 40 Small capacity CO ₂ Incubator, 230V, 50/60 Hz



© 2010 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries.

North America: USA/Canada +1 866 984 3766 (866-9-THERMO)

www.thermoscientific.com/incubators

Europe: Austria +43 1 801 40 0, Belgium +32 2 482 30 30, France +33 2 2803 2180, Germany national toll free 08001-536 376, Germany international +49 6184 90 6940, Italy +39 02 02 95059 434-254-375, Netherlands +31 76 571 4440, Nordic/Baltic countries +358 9 329 100, Russia/CIS +7 (812) 703 42 15, Spain/Portugal +34 93 223 09 18, Switzerland +41 44 454 12 12, UK/Ireland +44 870 609 9203

Asia: China +86 21 6865 4588 or +86 10 8419 3588, India toll free 1800 22 8374, India +91 22 6716 2200, Japan +81 45 453 9220, Other Asian countries +852 2885 4613 **Countries not listed:** +49 6184 90 6940 or +33 2 2803 2180