

Ultracentrifuge Rotors

Thermo Scientific Fiberlite F65L-6x13.5 Rotor

F65L-6x13.5



ROTOR SPECIFICATIONS

Type	Fixed Angle
Material	Carbon Fiber Composite
Capacity	6 x 13.5 mL
Maximum Speed	65,000 rpm
Tube Dimensions Ø x L	16 x 76 mm
Angle	30°
K Factor	70
Net Weight	4.31 kg/9.5 lb
Warranty	15 years

	RCF	RADIUS
Maximum	324,000 xg	68.4 mm
Minimum	99,372 xg	21 mm

ROTOR AND ACCESSORIES ORDERING INFORMATION

		THERMO SCIENTIFIC SORVALL WX Series/Discovery Series	BECKMAN™ Floor Model L8-Series and previous
Fiberlite F65L-6x13.5	Cat. No.	096-067135	096-067135
	Max. Speed	65,000 rpm	65,000 rpm
	Max. RCF	324,000 xg	324,000 xg

Fiberlite F65L-6x13.5 Rotor Package includes:

Description	Cat. No.
Rotor Lid Assembly	099-067135
O-ring Refresher Kit (Includes 5 O-rings, 3 Vacuum Grease Packets (65937), 3 Lubricant Packets (61556) and 4 oz FiberClene Lab Cleaner (017-0024))	021-067135
Lid Fastening Tool	014-2455
PC Bottle Assembly (qty. 3 pkg. of 2)	010-0128
Tube Removal Tool	097-1133
Rotor Stand	20048993
Vacuum Grease	65937
Lubricant	61556
Ultraspeed Rotor Log Book	52384
Extra Overspeed Decal (65,000 rpm)	014-2466

LABWARE ORDERING INFORMATION

Tube Vol. (mL)	Vol. Fill (mL)	Description	Cat. No.	Qty. /Pkg.	Tubes /Rotor	Max. Speed (rpm)	Max. RCF (xg)	Dims ØxL (mm)	REQUIRED TUBE CLOSURES, ADAPTERS, TOOLS & ACCESSORIES			
									Type	Cat. No.	Qty. /Pkg.	Description
13.5	10.4	PC Bottle Assembly	010-1382	2	6	65,000	324,000	16 x 76	Closure	Included	2	Aluminum Sealing Assembly
									Tool	097-1133	1	Removal Tool
13.5	10.4	PC Bottle Only	010-1380	6	6	65,000	324,000	16 x 76	—	—	—	—
13.5	8.9	OptiSeal™ Polyallomer Tube	010-1262	56	6	65,000	324,000	16 x 60	Closure	010-1355	2	Spacer
									Tool	017-0045	1	Spacer Removal Tool
									Tool	017-0044	1	Tube Extraction Tool

© 2014 Thermo Fisher Scientific Inc. All rights reserved. Beckman and OptiSeal are registered trademarks of Beckman Coulter Inc. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

For more information, please visit
www.thermoscientific.com/rotors.

DSCFGULTRAF65L6X13.5 1014

Thermo
 SCIENTIFIC

Part of Thermo Fisher Scientific